



CERTIFICATION

AOAC Research Institute *Performance Tested Methods*SM

Certificate No.
022003

The AOAC Research Institute hereby certifies the method known as

BACGene STEC TOP7 Workflow

manufactured by

Gold Standard Diagnostics Freiburg GmbH
Engesserstraße 4
D-79108 Freiburg im Breisgau
Germany

Gold Standard Diagnostics Budapest Kft.
Fóti út 56. Building A
HU-1047, Budapest
Hungary

This method has been evaluated and certified according to the policies and procedures of the AOAC *Performance Tested Methods*SM Program. This certificate indicates an AOAC Research Institute Certification Mark License Agreement has been executed which authorizes the manufacturer to display the AOAC Research Institute *Performance Tested Methods*SM certification mark on the above-mentioned method for the period below. Renewal may be granted by the Expiration Date under the rules stated in the licensing agreement.

A handwritten signature in black ink, appearing to read "Bradley A. Stawick".

Bradley A. Stawick, AOAC Research Institute Senior Director

Issue Date	December 19, 2025
Expiration Date	December 31, 2026

AUTHORS ORIGINAL VALIDATION: Nicole Klass, Benjamin Bastin, M. Joseph Benzinger, Jr., James Agin, Laura Bleichner, and Christoph Bahrtd MODIFICATION NOVEMBER 2020: Sándor Varga and Dr. Nikolett Horváth MODIFICATION JANUARY 2021: Laura Bleichner and Vera Bleicher MODIFICATION JANUARY 2024: Laura Bleichner and Dr. Stefan Kohl	SUBMITTING COMPANY Gold Standard Diagnostics Freiburg GmbH, formerly, Eurofins GeneScan Technologies GmbH Engesserstraße 4 D-79108 Freiburg im Breisgau Germany
--	---

METHOD NAME BACGene STEC TOP7 Workflow	CATALOG NUMBERS 1. BACGene <i>E. coli</i> O157:H7: 5123222701, 5123222711 2. BACGene Mplex STEC Screen: 5123223001, 5123223011 3. BACGene Mplex STEC SEROtype1: 5123223101, 5123223111 4. BACGene Mplex STEC SEROtype2: 5123223201, 5123223211
--	--

INDEPENDENT LABORATORY Q Laboratories 1930 Radcliff Drive Cincinnati, OH 45204 USA	APPLICABILITY OF METHOD Target organisms – STEC Top 7 (<i>Escherichia coli</i> O157:H7 and non-O157 STEC <i>E. coli</i> O26, O45, O103, O111, O121, O145). Matrixes – Fresh raw ground beef (375 g test portion), fresh raw beef trim (375 g test portion), beef carcass sampling cloth (4" x 4" cloth), Romaine lettuce (375 g test portion), all-purpose flour (25 g test portion), stainless steel (4" x 4" surface area, 18-gauge 300 series, brush finish, NSF-certified), unglazed ceramic tile (1" x 1" surface area). Performance claims – The BACGene STEC workflow is comparable to the U. S. Department of Agriculture Food Safety and Inspection Service (FSIS) <i>Microbiology Laboratory Guidebook</i> (MLG), 5C.00, <i>Detection, Isolation And Identification of Top Seven Shiga Toxin-Producing <i>Escherichia coli</i> (STECs) from Meat Products and Carcass and Environmental Sponges</i> (2) (fresh raw ground beef (73% lean), fresh raw beef trim, beef carcass sampling cloth, stainless steel environmental sponges and ceramic tile environmental swabs) and the Food and Drug Administration (FDA) <i>Bacteriological Analytical Manual</i> (BAM), Chapter 4A, <i>Diarrheagenic <i>Escherichia coli</i></i> (3) reference methods (Romaine lettuce and all-purpose flour).
--	--

ORIGINAL CERTIFICATION DATE February 08, 2020	CERTIFICATION RENEWAL RECORD Renewed through December 2026.
METHOD MODIFICATION RECORD 1. November 2020 Level 2 2. January 2021 Level 2 3. December 2022 Level 1 4. January 2024 Level 2 5. February 2025 Level 1	SUMMARY OF MODIFICATION 1. Manufacturing location addition of Eurofins Technologies Hungary Kft. 2. Detergent replacement required by European Chemicals Agency and edits. 3. Company rebranded to Gold Standard Diagnostics Freiburg GmbH. 4. Addition of PURE software. 5. Editorial changes.
Under this AOAC <i>Performance Tested Methods</i> SM License Number, 022002 this method is distributed by: NONE	Under this AOAC <i>Performance Tested Methods</i> SM License Number, 022002 this method is distributed as: NONE

PRINCIPLE OF THE METHOD (1)

The BACGene STEC TOP7 workflow utilizes the BACGene Mplex STEC Screen, BACGENE *E. coli* O157:H7, BACGene Mplex SEROtype 1 and BACGene Mplex SEROtype 2 test kits for the detection of the "Top 7" Shiga-Toxin Producing *E. coli* (O157, O26, O45, O103, O111, O121, and O145) containing *stx1* and/or *stx2* and *eae*. The BACGene technology is based on the detection of amplified fragments using FAM™ fluorescence-labeled hybridization probe quenched by non-fluorescent Tide Quencher™ 2 or by non-fluorescent Quencher MGBEQ; R6G™ fluorescence-labeled hybridization probe quenched by non-fluorescent Tide Quencher 2; or ROX™ fluorescence-labeled hybridization probe quenched by non-fluorescent Tide Quencher 3.

The BACGene Mplex STEC workflow utilizes Uracil-DNA Glycosylase (UNG) where Deoxyuridine triphosphate (dUTP) substitutes for deoxythymidine triphosphate (dTTP) in the PCR master mix. This substitution causes newly synthesized amplicons to contain deoxyuridine instead of deoxythymidine. The heat labile UNG cleaves hydrolytically at the N-glycosylic bond of uracil and deoxyribose where it is inactivated by the initial denaturation of the PCR. A PCR preceded by this step exclusively amplifies from native, i.e. thymine containing DNA (e.g. bacterial DNA). Thus, reducing the risk of false positive signals due to carry over contamination from previous analyses.

The BACGene STEC workflow consists of a screening step with the BACGene Mplex STEC Screen test kit followed by a screening step with the BACGene *E. coli* O157:H7, BACGene Mplex SEROtype 1 and BACGene Mplex SEROtype 2 test kits. A test sample is deemed a presumptive positive for the identified *E. coli* serotype if the sample receives a positive result from both the BACGene Mplex STEC Screen kit and the BACGene *E. coli* O157:H7, BACGene Mplex SEROtype 1 or BACGene Mplex SEROtype 2 test kits. This procedure allows for the detection of the virulence genes first and then the detection of a specific serotype, thus allowing the user to detect any STEC contamination present.

The BACGene STEC TOP7 workflow was validated for use with the Agilent AriaMx™ and/or Bio-Rad CFX96 Touch™ Deep Well PCR platforms. All run files were then downloaded and analyzed by either the BACGene Evaluation Sheet or the FastFinder software. Additionally, there is an option of conducting a sample pre-treatment using PRERaser BACGene to eliminate free DNA from enrichments prior to the lysis procedure.

All test portions are enriched in pre-warmed modified Tryptone Soya Broth (mTSB) for 10-24 h at 41.5±1°C. After enrichment, aliquots of 30 µL is sampled and thermally and enzymatically lysed to release the DNA and 5 µL of the lysate is then analyzed by real-time PCR using Agilent AriaMx™ and/or Bio-Rad CFX96 Touch™ Deep Well PCR platforms.

DISCUSSION OF THE VALIDATION STUDY (1)

BACGene Mplex STEC Workflow

The BACGene Mplex STEC workflow successfully recovered the Top 7 STECs, *E. coli* O157:H7, O26, O45, O103, O111, O121, and O145, from fresh raw ground beef, fresh raw beef trim, beef carcass sampling cloth, romaine lettuce, all-purpose flour, stainless steel environmental sponges and ceramic tile environmental swabs after 10 and 24 hours of enrichment. Using POD analysis, no statistically significant differences were observed between the number of positive samples detected by the BACGene Mplex STEC workflow and the reference method for all samples tested for both timepoints regardless of PCR platform or the use of the PRERaser BACGene treatment.

Additionally, the alternative confirmation conducted according to the Abraxis IMS beads for BACGene *E. coli* STEC Solution workflow successfully recovered each Top 7 STEC from each matrix. The BACGene Mplex STEC Screen, Mplex SEROtype1 and Mplex SEROtype2 test kits accurately detected the target analytes and correctly exclude all non-target organisms.

The BACGene STEC TOP7 workflow is quick to perform, on the Agilent AriaMx™ and the Bio-Rad CFX96 Touch Deep Well, with or without the PRERaser BACGene pre-treatment. The Abraxis IMS beads for BACGene *E. coli* STEC Solution workflow utilizes a small footprint and less consumables than the reference method procedures. Confirmation results can be provided in less than 2 hours post-incubation for 94 sample replicates when using either workflow. By being capable of being ran on multiple PCR platforms and the optional use of a pre-treatment, the BACGene STEC TOP7 workflow is customizable to the end user.

NOTE: Data generated for the Inclusivity, Exclusivity, and Matrix studies were generated by Q Laboratories.

Data Tables for Inclusivity and Exclusivity (1)

Table 54. BACGene *E. coli* O157:H7 Kit, Inclusivity Results (1)

No.	Species	Serotype	Source	Origin	Result	
					Agilent AriaMx	Bio-Rad CFX96 Touch Deep Well
1	<i>E. coli</i>	O157:H7	MSU ¹ TW00116	Human	+	+
2	<i>E. coli</i>	O157:H7	MSU TW00975	Human	+	+
3	<i>E. coli</i>	O157:H7	MSU TW02302	Hamburger	+	+
4	<i>E. coli</i>	O157:H7	MSU TW04863	Human	+	+
5	<i>E. coli</i>	O157:H7	MSU TW05356	Human	+	+
6	<i>E. coli</i>	O157:H7	MSU TW07587	Human	+	+
7	<i>E. coli</i>	O157:H7	ATCC ² BAA-460	Human Feces	+	+
8	<i>E. coli</i>	O157:H7	ATCC 35150	Human Feces	+	+
9	<i>E. coli</i>	O157:H7	ATCC 43889	Human Feces	+	+
10	<i>E. coli</i>	O157:H7	ATCC 43895	Hamburger	+	+
11	<i>E. coli</i>	O157:H7	ATCC 51657	Clinical Isolate	+	+
12	<i>E. coli</i>	O157:H7	ATCC 51659	Clinical Isolate	+	+
13	<i>E. coli</i>	O157:H7	ATCC 700599	Salami	+	+
14	<i>E. coli</i>	O157:H7	ATCC 700927	Not Available	+	+
15	<i>E. coli</i>	O157:H7	MSU DEC3A	Human	+	+
16	<i>E. coli</i>	O157:H7	MSU DEC3B	Human	+	+
17	<i>E. coli</i>	O157:H7	MSU DEC3C	Human	+	+
18	<i>E. coli</i>	O157:H7	MSU DEC3D	Human	+	+
19	<i>E. coli</i>	O157:H7	MSU DEC3E	Human	+	+
20	<i>E. coli</i>	O157:H7	MSU DEC4B	Human	+	+
21	<i>E. coli</i>	O157:H7	MSU DEC4D	Cow, Calf	+	+
22	<i>E. coli</i>	O157:H7	MSU DEC4E	Human	+	+
23	<i>E. coli</i>	O157:H7	QL ³ 2-202	Meat	+	+
24	<i>E. coli</i>	O157:H7	QL 2-205	Meat	+	+
25	<i>E. coli</i>	O157:H7	QL 2-206	Meat	+	+
26	<i>E. coli</i>	O157:H7	QL 2-207	Meat	+	+

27	<i>E. coli</i>	O157:H7	QL 2-370	Meat	+	+
28	<i>E. coli</i>	O157:H7	QL 2-701	Beef	+	+
29	<i>E. coli</i>	O157:H7	QL 2-704	Beef	+	+
30	<i>E. coli</i>	O157:H7	QL 2-705	Beef	+	+
31	<i>E. coli</i>	O157:H7	QL 2-706	Beef	+	+
32	<i>E. coli</i>	O157:H7	QL 2-707	Beef	+	+
33	<i>E. coli</i>	O157:H7	USDA ⁴ 465-97	Beef	+	+
34	<i>E. coli</i>	O157:H7	QL 2-709	Beef	+	+
35	<i>E. coli</i>	O157:H7	QL 2-710	Beef	+	+
36	<i>E. coli</i>	O157:H7	QL 14077.2	Meat	+	+
37	<i>E. coli</i>	O157:H7	QL 14077.3	Salami	+	+
38	<i>E. coli</i>	O157:H7	QL 14077.5	Hamburger	+	+
39	<i>E. coli</i>	O157:H7	QL 14077.6	Food	+	+

Table 54. BACGene *E. coli* O157:H7 Kit, Inclusivity Results (Continued) (1)

No.	Species	Serotype	Source	Origin	Result	
					Agilent AriaMx	Bio-Rad CFX96 Touch Deep Well
40	<i>E. coli</i>	O157:H7	QL 14077.7	Ground Beef	+	+
41	<i>E. coli</i>	O157:H7	QL 14077.8	Ground Beef	+	+
42	<i>E. coli</i>	O157:H7	ATCC 700376	Human Feces	+	+
43	<i>E. coli</i>	O157:H7	BEI NR ⁵ -4356	Not Available	+	+
44	<i>E. coli</i>	O157:H7	ATCC 43890	Human Feces	+	+
45	<i>E. coli</i>	O157:H7	QL 12298-1	Not Available	+	+
46	<i>E. coli</i>	O157:H7	QL 12298-2	Deli Ham	+	+
47	<i>E. coli</i>	O157:H7	QL 14077.1	Meat	+	+
48	<i>E. coli</i>	O157:H7	QL 2-214	Meat	+	+
49	<i>E. coli</i>	O157:H7	BEI NR-12	Not Available	+	+
50	<i>E. coli</i>	O157:H7	QL 2-203	Meat	+	+

1. MSU – Michigan State University Culture Collection, 2. ATCC – American Type Culture Collection, 3. QL – Q Laboratories Culture Collection, 4. USDA – United States Department of Agriculture Culture Collection, 5. BEI-NR - BEI Resources Culture Collection

Table 55. BACGene *E. coli* O157:H7 Kit, Exclusivity Results (1)

No.	Species	Source	Origin	Result	
				Agilent AriaMx	Bio-Rad CFX96 Touch Deep Well
1	<i>Aeromonas hydrophila</i>	ATCC ¹ 49140	Clinical Isolate	-	-
2	<i>Citrobacter braakii</i>	ATCC 43162	Clinical Isolate	-	-
3	<i>Citrobacter farmeri</i>	ATCC 51633	Human Feces	-	-
4	<i>Cronobacter sakazakii</i>	QL ² 17031.4	Infant Formula	-	-
5	<i>E. coli</i> O103	MSU ³ TW08101	Human	-	-
6	<i>E. coli</i> O103	PSU ⁴ 9.0036	Not Available	-	-
7	<i>E. coli</i> O111	MSU DEC8D	Human Infant	-	-
8	<i>E. coli</i> O121	PSU 5.0959	Not Available	-	-
9	<i>E. coli</i> O121	PSU 7.1686	Not Available	-	-
10	<i>E. coli</i> O145	MSU TW07596	Human	-	-
11	<i>E. coli</i> O26	ATCC BAA-1653	Stool, Associated with Ice Cream	-	-
12	<i>E. coli</i> O26	MSU TW07814	Human	-	-
13	<i>E. coli</i> O45	MSU TW14003	Human	-	-

14	<i>E. coli</i> O45	MSU TW07947	Human	-	-
15	<i>Edwardsiella tarda</i>	ATCC 15947	Human Feces	-	-
16	<i>Enterobacter aerogenes</i>	ATCC 13048	Sputum	-	-
17	<i>Escherichia fergusonii</i>	ATCC 35469	Human Feces	-	-
18	<i>Escherichia hermanii</i>	ATCC 33650	Mouse Brain	-	-
19	<i>Escherichia vulneris</i>	ATCC 29943	Clinical Isolate, Human Wound	-	-
20	<i>Haemophilus influenzae</i>	ATCC 19418	Not Available	-	-
21	<i>Hafnia alvei</i>	ATCC 51815	Milk	-	-
22	<i>Klebsiella pneumoniae</i>	ATCC 4352	Dairy Product	-	-
23	<i>Morganella morganii</i>	ATCC 25829	Human	-	-
24	<i>Mycobacterium smegmatis</i>	ATCC 19420	Not Available	-	-
25	<i>Pantoea agglomerans</i>	ATCC 19552	Sewage	-	-
26	<i>Providencia rettgeri</i>	ATCC 14505	Not Available	-	-
27	<i>Pseudomonas aeruginosa</i>	ATCC 9027	Clinical Isolate, Outer Ear Infection	-	-
28	<i>Serratia marcescens</i>	ATCC 13880	Human	-	-
29	<i>Shigella boydii</i>	ATCC 9290	Human Feces	-	-
30	<i>Vibrio vulnificus</i>	ATCC 49140	Clinical Isolate	-	-

ATCC – American Type Culture Collection, 2. QL – Q Laboratories Culture Collection, 3. MSU – Michigan State University Culture Collection, 4. PSU - Pennsylvania State University Culture Collection

Table 56. BACGene Mplex STEC Screen Kit, Inclusivity Results (1)

No.	Species and Serotype	Source	Known Virulence Factors			Origin	Result					
							Agilent AriaMx			Bio-Rad CFX96 Touch Deep Well		
			stx1	stx2	eae		stx1	stx2	eae	stx1	stx2	eae
1	<i>E. coli</i> O103	MSU ¹ TW08101	+	-	+	Human	+	-	+	+	-	+
2	<i>E. coli</i> O103	MSU TW04162	+	-	+	Human	+	-	+	+	-	+
3	<i>E. coli</i> O103	MSU TW05997	+	-	+	Human	+	-	+	+	-	+
4	<i>E. coli</i> O103	MSU TW07697	+	-	+	Human	+	-	+	+	-	+
5	<i>E. coli</i> O103	MSU TW07971	+	-	+	Human	+	-	+	+	-	+
6	<i>E. coli</i> O103	MSU TW11239	+	-	+	Human	+	+ ^a / _{-b}	+	+	-	+
7	<i>E. coli</i> O103	NCTC ² 8196	/ ^c	/ ^c	+	Not Available	+	-	+	+	-	+
8	<i>E. coli</i> O103	PSU ³ 9.0036	/ ^c	/ ^c	+	Not Available	+	+	+	+	+	+
9	<i>E. coli</i> O111	MSU DEC 8D	+	+	+	Human Infant	+	+	+	+	+	+
10	<i>E. coli</i> O111	MSU DEC15A	/ ^c	/ ^c	+	Human	+	+	+	+	+	+
11	<i>E. coli</i> O111	MSU TW00186	+	+	+	Human	+	+	+	+	+	+
12	<i>E. coli</i> O111	MSU TW01387	+	-	+	Human	+	-	+	+	-	+
13	<i>E. coli</i> O111	MSU TW05614	+	+	+	Human	+	+	+	+	+	+
14	<i>E. coli</i> O111	MSU TW06315	+	+	+	Human	+	+	+	+	+	+
15	<i>E. coli</i> O111	MSU TW07502	+	-	+	Human	+	-	+	+	-	+
16	<i>E. coli</i> O111	MSU TW14960	+	-	+	Human	+	-	+	+	-	+
17	<i>E. coli</i> O111:H8	QL ⁴ 12289-3A	+	-	+	Not Available	+	-	+	+	-	+
18	<i>E. coli</i> O121	MSU MT#2	-	+	+	Human	+	+	+	+	+	+
19	<i>E. coli</i> O121	MSU TW07614	-	+	+	Human	+	+	+	+	+	+
20	<i>E. coli</i> O121	MSU TW07931	-	+	+	Human	-	+	+	+	+	+
21	<i>E. coli</i> O121	MSU TW08039	-	+	+	Human	-	+	+	-	+	+
22	<i>E. coli</i> O121	PSU 10.0709	/ ^c	/ ^c	+	Not Available	+	+ ^a / _{-b}	+	+ ^a / _{-b}	+	+
23	<i>E. coli</i> O121	PSU 5.0959	-	+	+	Not Available	-	+	+	-	+	+
24	<i>E. coli</i> O121	PSU 7.1686	/ ^c	/ ^c	+	Not Available	+	+	+	+	+	+
25	<i>E. coli</i> O121	PSU 7.1709	/ ^c	/ ^c	+	Not Available	+	+	+	+	+	+
26	<i>E. coli</i> O121	PSU 7.1732	/ ^c	/ ^c	+	Not Available	+	+	+	+	+	+
27	<i>E. coli</i> O145	MSU TW01664	+	-	+	Human	+	+	+	+	+	+
28	<i>E. coli</i> O145	MSU TW07596	+	-	+	Human	+	-	+	+	-	+
29	<i>E. coli</i> O145	MSU TW09153	+	-	+	Human	+	-	+	+	-	+
30	<i>E. coli</i> O145	MSU TW09356	+	-	+	Human	+ ^a / _{-b}	+	+	+	-	+
31	<i>E. coli</i> O145	PSU 10.0707	/ ^c	/ ^c	+	Not Available	+	+	+	+	+	+
32	<i>E. coli</i> O145	PSU 10.0708	/ ^c	/ ^c	+	Not Available	+	-	+	+	+	+
33	<i>E. coli</i> O145	PSU 7.1711	/ ^c	/ ^c	+	Not Available	-	+	+	-	+	+
34	<i>E. coli</i> O26	ATCC ⁵ BAA-1653	+	+	+	Human Feces	+	-	+	+	-	+
35	<i>E. coli</i> O26	MSU DEC 10B	+	-	+	Human	+	-	+	+	-	+
36	<i>E. coli</i> O26	MSU DEC 9F	+	-	+	Human	+	-	+	+	-	+
37	<i>E. coli</i> O26	MSU TW07862	+	-	+	Calf, Cow	+	-	+	+	-	+
38	<i>E. coli</i> O26	MSU TW00971	+	-	+	Human Feces	+	-	+	+	-	+
39	<i>E. coli</i> O26	MSU TW02295	+	-	+	Human Infant	+	-	+	+	-	+
40	<i>E. coli</i> O26	MSU TW04284	+	-	+	Human	+	-	+	+	-	+
41	<i>E. coli</i> O26	MSU TW07814	+	+	+	Human	+	+	+	+	+	+

42	<i>E. coli</i> O26	MSU TW08031	+	-	+	Human	+	-	+	+	-	+
43	<i>E. coli</i> O45	MSU DEC 11C	/ ^c	/ ^c	+	Human	+	-	+	+	-	+
44	<i>E. coli</i> O45	MSU TW07947	+	-	+	Human	+	-	+	+	-	+
45	<i>E. coli</i> O45	MSU TW09183	+	-	+	Human	+	-	+	+	-	+
46	<i>E. coli</i> O45	MSU TW10121	+	-	+	Human	+	-	+	+	-	+
47	<i>E. coli</i> O45	MSU TW14003	+	-	+	Human	+	-	+	+	-	+
48	<i>E. coli</i> O157:H7	BEI NR-4356	/ ^c	/ ^c	+	Not Available	+ ^a / ^{-b}	+	+	-	+	+
49	<i>E. coli</i> O157:H7	MSU TW00116	-	+	+	Not Available	-	+	+	+ ^a / ^{-b}	+	+
50	<i>E. coli</i> O157:H7 (NM)	ATCC 700376	+	-	+	Not Available	+	-	+	+	-	+

1. MSU – Michigan State University Culture Collection, 2. NCTC – National Collection of Type Cultures, 3. PSU - Pennsylvania State University Culture Collection, 4.

QL – Q Laboratories Culture Collection, 5. ATCC – American Type Culture Collection

^aResults are positive when eae as well as one of the shigatoxins (stx1 or stx2) is positive

Table 57. BACGene Mplex STEC Screen Kit, Exclusivity Results (1)

No.	Species	Source	Origin	Result	
				Agilent AriaMx	Bio-Rad CFX96 Touch Deep Well
1	<i>Aeromonas hydrophila</i>	ATCC ¹ 49140	Clinical Isolate	-	-
2	<i>Alcaligenes faecalis</i>	ATCC 8750	Not Available	-	-
3	<i>Citrobacter braakii</i>	ATCC 43162	Clinical Isolate	-	-
4	<i>Citrobacter farmeri</i>	ATCC 51633	Human Feces	- ^a	-
5	<i>Cronobacter sakazakii</i>	QL ² 17031.4	Infant Formula	-	-
6	<i>Escherichia coli</i>	ATCC 8739	Feces	-	-
7	<i>Edwardsiella tarda</i>	ATCC 15947	Human Feces	-	-
8	<i>Enterobacter aerogenes</i>	ATCC 35029	Not Available	-	-
9	<i>Escherichia blattae</i>	ATCC 29907	Insect	-	-
10	<i>Escherichia fergusonii</i>	ATCC 35469	Human Feces	-	-
11	<i>Escherichia hermannii</i>	ATCC 33650	Mouse Brain	-	-
12	<i>Escherichia vulneris</i>	ATCC 29943	Clinical Isolate, Human Wound	-	-
13	<i>Haemophilus influenzae</i>	ATCC 19418	Not Available	-	-
14	<i>Hafnia alvei</i>	ATCC 51815	Milk	-	-
15	<i>Klebsiella pneumoniae</i>	ATCC 13440	Not Available	-	-
16	<i>KLuyvera intermedia</i>	ATCC 33110	Surface Water	-	-
17	<i>Morganella morganii</i>	ATCC 25829	Human	-	-
18	<i>Mycobacterium smegmatis</i>	ATCC 19420	Not Available	-	-
19	<i>Pantoea agglomerans</i>	ATCC 19552	Sewage	-	-
20	<i>Proteus mirabilis</i>	ATCC 7002	Urine	-	-
21	<i>Providencia rettgeri</i>	ATCC 14505	Not Available	-	-
22	<i>Pseudomonas aeruginosa</i>	ATCC 9027	Clinical Isolate, Outer Ear Infection	-	-
23	<i>Rahnella aquatilis</i>	ATCC 55046	Soil	-	-
24	<i>Salmonella Typhimurium</i>	ATCC 14028	Animal Tissue	-	-
25	<i>Serratia marcescens</i>	ATCC 8100	Not Available	-	-
26	<i>Shigella boydii</i>	ATCC 9207	Pork Liver	-	-
27	<i>Siccibacter turicensis</i>	ATCC 54945	Not Available	-	-
28	<i>Sphingomonas paucimobilis</i>	ATCC 29837	Hospital Respirator	- ^b	-
29	<i>Vibrio vulnificus</i>	QL 021111A	Seafood Product	-	-
30	<i>Yersinia enterocolitica</i>	ATCC 49397	Clinical Isolate	-	-

¹. ATCC – American Type Culture Collection, ². QL – Q Laboratories Culture Collection

^aIsolate produced a positive result for “stx2” when analyzed by the FastFinder Software. ^bIsolate produced a positive result for “eae” when analyzed by the FastFinder Software.

Table 58. BACGene Mplex SEROtype1 Kit, Inclusivity Results (1)

No.	Species	Serotype	Source	Origin	Result	
					Agilent AriaMx	Bio-Rad CFX96 Touch Deep Well
1	<i>E. coli</i>	O26	ATCC ¹ BAA-1653	Human Feces	+	+
2	<i>E. coli</i>	O26	MSU ² TW07862	Calf, Cow	+	+
3	<i>E. coli</i>	O26	MSU TW00971	Human Feces	+	+
4	<i>E. coli</i>	O26	MSU TW04284	Human	+	+
5	<i>E. coli</i>	O26	MSU TW08031	Human	+	+
6	<i>E. coli</i>	O26	MSU TW02295	Human Infant	+	+
7	<i>E. coli</i>	O26	MSU TW07814	Human	+	+
8	<i>E. coli</i>	O26	MSU TW05992	Human	+	+
9	<i>E. coli</i>	O26	MSU DEC 9A	Human	+	+
10	<i>E. coli</i>	O26	MSU DEC 9F	Human	+	+
11	<i>E. coli</i>	O26	MSU DEC 9F	Human	+	+
12	<i>E. coli</i>	O26	MSU DEC 10B	Human	+	+
13	<i>E. coli</i>	O26	MSU DEC 10C	Human	+	+
14	<i>E. coli</i>	O26:H11	MSU DEC 10E	Calf, Cow	+	+
15	<i>E. coli</i>	O103	MSU TW09101	Human	+	+
16	<i>E. coli</i>	O103	MSU TW07971	Human	+	+
17	<i>E. coli</i>	O103	MSU TW11239	Human	+	+
18	<i>E. coli</i>	O103	MSU TW07697	Human	+	+
19	<i>E. coli</i>	O103	MSU TW05997	Human	+	+
20	<i>E. coli</i>	O103	MSU TW04162	Human	+	+
21	<i>E. coli</i>	O103	MSU TW08101	Human	+	+
22	<i>E. coli</i>	O103	NCTC ³ 8196	Not Available	+	+
23	<i>E. coli</i>	O103	PSU ⁴ 5.1658	Not Available	+	+
24	<i>E. coli</i>	O103	PSU 7.1691	Not Available	+	+
25	<i>E. coli</i>	O103	PSU 9.0036	Not Available	+	+
26	<i>E. coli</i>	O103	QL ⁵ 15071-2	Meat Powder	+	+
27	<i>E. coli</i>	O111	MSU TW07926	Human	+	+
28	<i>E. coli</i>	O111	MSU TW14960	Human	+	+
29	<i>E. coli</i>	O111	MSU TW06296	Human	+	+
30	<i>E. coli</i>	O111	MSU TW06315	Human	+	+
31	<i>E. coli</i>	O111	MSU TW05614	Human	+	+
32	<i>E. coli</i>	O111	MSU TW00186	Human	+	+
33	<i>E. coli</i>	O111	MSU TW01387	Human	+	+
34	<i>E. coli</i>	O111	MSU TW07502	Human	+	+
35	<i>E. coli</i>	O111	MSU DEC 8D	Human Infant	+	+
36	<i>E. coli</i>	O111	MSU DEC 12A	Human Infant	+	+
37	<i>E. coli</i>	O111	MSU DEC15A	Human	+	+
38	<i>E. coli</i>	O111:H12	MSU DEC6A	Human Infant	+	+
39	<i>E. coli</i>	O111:H8	MSU DEC6C	Human	+	+
40	<i>E. coli</i>	O111:H8	QL 12289-3A	Not Available	+	+
41	<i>E. coli</i>	O26	FDA ⁶ 1143232-005	Food	+	+
42	<i>E. coli</i>	O26	FDA 1143052-008	Food	+	+
43	<i>E. coli</i>	O111	FDA 1146086-007	Food	+	+
44	<i>E. coli</i>	O103	FDA 1145997-002	Food	+	+
45	<i>E. coli</i>	O26	FDA AM00324	Clinical Isolate	+	+
46	<i>E. coli</i>	O26:K:H1	Eurofins ⁷	Not Available	+	+
47	<i>E. coli</i>	O103:H2	Eurofins	Not Available	+	+
48	<i>E. coli</i>	O111K:H-	Eurofins	Not Available	+	+
49	<i>E. coli</i>	O111:H8	Eurofins	Not Available	+	+
50	<i>E. coli</i>	O26:H11	Eurofins	Not Available	+	+
51	<i>E. coli</i>	O111:H-	Eurofins	Not Available	+	+

1. ATCC – American Type Culture Collection, 2. MSU – Michigan State University Culture Collection, 3. NCTC – National Collection of Type Cultures, 4. PSU – Pennsylvania State University Culture Collection, 5. QL – Q Laboratories Culture Collection, 6. FDA – U.S. Food and Drug Administration Culture Collection, 7. Eurofins – Eurofins Culture Collection (Isolates with this source identification were analyzed by Eurofins GeneScan Technologies, GmbH in Germany)

Table 59. BACGene Mplex SEROtype1 Kit, Exclusivity Results (1)

No.	Species	Source	Origin	Result	
				Agilent AriaMx	Bio-Rad CFX96 Touch Deep Well
1	<i>E. coli</i> O45	MSU ¹ TW14003	Human	-	-
2	<i>E. coli</i> O45	MSU TW07947	Human	-	-
3	<i>E. coli</i> O45	MSU DEC 11C	Human	-	-
4	<i>E. coli</i> O121	PSU ² 5.0959	Not Available	-	-

5	<i>E. coli</i> O121	MSU TW07931	Human	-	-
6	<i>E. coli</i> O121	MSU TW08039	Human	-	-
7	<i>E. coli</i> O145	PSU 7.1711	Not Available	-	-
8	<i>E. coli</i> O145	PSU 10.0707	Not Available	-	-
9	<i>E. coli</i> O145	MSU TW07596	Human	-	-
10	<i>E. coli</i> O145	MSU TW09153	Human	-	-
11	<i>Aeromonas hydrophila</i>	ATCC ³ 49140	Clinical Isolate	-	-
12	<i>Alcaligenes faecalis</i>	ATCC 8750	Not Available	-	-
13	<i>Cronobacter sakazakii</i>	QL ⁴ 17031.4	Infant Formula	-	-
14	<i>Edwardsiella tarda</i>	ATCC 15947	Human Feces	-	-
15	<i>Escherichia blattae</i>	ATCC 29907	Insect	-	-
16	<i>Escherichia hermarii</i>	ATCC 33650	Mouse Brain	-	-
17	<i>Escherichia vulneris</i>	ATCC 29943	Clinical Isolate, Human Wound	-	-
18	<i>Haemophilus influenzae</i>	ATCC 19418	Not Available	-	-
19	<i>Hafnia alvei</i>	ATCC 51815	Milk	-	-
20	<i>Klebsiella pneumoniae</i>	ATCC 4352	Not Available	-	-
21	<i>Kluyvera intermedia</i>	ATCC 33110	Surface Water	-	-
22	<i>Mycobacterium smegmatis</i>	ATCC 19420	Not Available	-	-
23	<i>Pantoea agglomerans</i>	ATCC 19552	Sewage	-	-
24	<i>Proteus mirabilis</i>	ATCC 7002	Urine	-	-
25	<i>Providencia rettgeri</i>	ATCC 14505	Not Available	-	-
26	<i>Serratia marcescens</i>	ATCC 13880	Not Available	-	-
27	<i>Shigella boydii</i>	ATCC 9290	Pork Liver	-	-
28	<i>Siccibacter turicensis</i>	CCUG ⁵ 54945	Not Available	-	-
29	<i>Sphingomonas paucimobilis</i>	ATCC 29837	Hospital Respirator	-	-
30	<i>Vibrio vulnificus</i>	QL 02111-1A	Seafood Product	-	-

1. MSU – Michigan State University Culture Collection, 2. PSU - Pennsylvania State University Culture Collection, 3. ATCC – American Type Culture Collection, 4. QL – Q Laboratories Culture Collection, 5. CCUG - Culture Collection University of Gothenburg

Table 60. BACGene Mplex SEROtype2 Kit, Inclusivity Results (1)

No.	Species	Serotype	Source	Origin	Result	
					Agilent AriaMx	Bio-Rad CFX96 Touch Deep Well
1	<i>E. coli</i>	O45	MSU1 TW10121	Human	+	+
2	<i>E. coli</i>	O45	MSU TW14003	Human	+	+
3	<i>E. coli</i>	O45	MSU TW07947	Human	+	+
4	<i>E. coli</i>	O45	MSU DEC 11C	Human	+	+
5	<i>E. coli</i>	O45	MSU TW09183	Human	+	+
6	<i>E. coli</i>	O45	PSU ² 11.1079	Not Available	+	+
7	<i>E. coli</i>	O45	PSU 1.2635	Not Available	+	+
8	<i>E. coli</i>	O45	PSU 6.3127	Not Available	+	+
9	<i>E. coli</i>	O121	MSU TW08039	Human	+	+
10	<i>E. coli</i>	O121	PSU 10.0709	Not Available	+	+
11	<i>E. coli</i>	O121	PSU 5.0959	Not Available	+	+
12	<i>E. coli</i>	O121	PSU 7.1686	Not Available	+	+
13	<i>E. coli</i>	O121	PSU 7.1709	Not Available	+	+
14	<i>E. coli</i>	O121	PSU 7.1732	Not Available	+	+
15	<i>E. coli</i>	O121	MSU TW07931	Human	+	+
16	<i>E. coli</i>	O121	MSU TW07614	Human	+	+
17	<i>E. coli</i>	O121	MSU MT#2	Human	+	+
18	<i>E. coli</i>	O145	PSU 7.1711	Not Available	+	+
19	<i>E. coli</i>	O145	PSU 10.0707	Not Available	+	+
20	<i>E. coli</i>	O145	PSU 10.0708	Not Available	+	+
21	<i>E. coli</i>	O145	MSU TW07596	Human	+	+
22	<i>E. coli</i>	O145	MSU TW09356	Human	+	+
23	<i>E. coli</i>	O145	MSU TW01664	Human	+	+
24	<i>E. coli</i>	O145	MSU TW09153	Human	+	+
25	<i>E. coli</i>	O121	MSU TW08023	Human	+	+
26	<i>E. coli</i>	O145	FDA ³ AM00376	Clinical Isolate	+	+
27	<i>E. coli</i>	O145	FDA AM00297	Clinical Isolate	+	+
28	<i>E. coli</i>	O45	FDA AM00274	Clinical Isolate	+	+
29	<i>E. coli</i>	O45	FDA AM00275	Clinical Isolate	+	+
30	<i>E. coli</i>	O45	FDA AM00276	Clinical Isolate	+	+
31	<i>E. coli</i>	O45	FDA AM00285	Clinical Isolate	+	+
32	<i>E. coli</i>	O45	FDA AM00291	Clinical Isolate	+	+
33	<i>E. coli</i>	O45	FDA AM00296	Clinical Isolate	+	+
34	<i>E. coli</i>	O45	FDA AM00304	Clinical Isolate	+	+
35	<i>E. coli</i>	O45	FDA AM00308	Clinical Isolate	+	+

36	<i>E. coli</i>	O45	FDA AM00311	Clinical Isolate	+	+
37	<i>E. coli</i>	O45	FDA AM00318	Clinical Isolate	+	+
38	<i>E. coli</i>	O45	FDA AM00319	Clinical Isolate	+	+
39	<i>E. coli</i>	O45	FDA AM00325	Clinical Isolate	+	+
40	<i>E. coli</i>	O45	FDA AM00328	Clinical Isolate	+	+
41	<i>E. coli</i>	O45	FDA AM00334	Clinical Isolate	+	+
42	<i>E. coli</i>	O45	FDA AM00335	Clinical Isolate	+	+
43	<i>E. coli</i>	O121	FDA AM00321	Clinical Isolate	+	+
44	<i>E. coli</i>	O45	FDA AM00345	Clinical Isolate	+	+
45	<i>E. coli</i>	O121	FDA AM00369	Clinical Isolate	+	+
46	<i>E. coli</i>	O45	FDA AM00347	Clinical Isolate	+	+
47	<i>E. coli</i>	O45	Eurofins ⁴	Not Available	+	+
48	<i>E. coli</i>	O45:H2	Eurofins	Not Available	+	+
49	<i>E. coli</i>	O121	Eurofins	Not Available	+	+
50	<i>E. coli</i>	O145	Eurofins	Not Available	+	+

1. MSU – Michigan State University Culture Collection, 2. PSU - Pennsylvania State University Culture Collection, 3. FDA – U.S. Food and Drug Administration Culture Collection, 4. Eurofins – Eurofins Culture Collection (Isolates with this source identification were analyzed by Eurofins GeneScan Technologies, GmbH in Germany)

Table 61. BACGene Mplex SEROtype2 Kit, Exclusivity Results (1)

No.	Species	Source	Origin	Result	
				Agilent AriaMx	Bio-Rad CFX96 Touch Deep Well
1	<i>E. coli</i> O26	ATCC ¹ BAA-1653	Human feces	-	-
2	<i>E. coli</i> O26	MSU ² TW07814	Calf, Cow	-	-
3	<i>E. coli</i> O26	MSU TW00971	Human Feces	-	-
4	<i>E. coli</i> O26	MSU TW04284	Human	-	-
5	<i>E. coli</i> O103	MSU TW09101	Human	-	-
6	<i>E. coli</i> O103	NCTC ³ 8196	Not Available	-	-
7	<i>E. coli</i> O103	PSU ⁴ 9.0036	Not Available	-	-
8	<i>E. coli</i> O111:H8	QL ⁵ 12289-3A	Not Available	-	-
9	<i>E. coli</i> O111	MSU TW01387	Human	-	-
10	<i>E. coli</i> O111	MSU TW07502	Human	-	-
11	<i>Aeromonas hydrophila</i>	ATCC 49140	Clinical Isolate	-	-
12	<i>Alcaligenes faecalis</i>	ATCC 8750	Not Available	-	-
13	<i>Cronobacter sakazakii</i>	QL 17031.4	Infant Formula	-	-
14	<i>Edwardsiella tarda</i>	ATCC 15947	Human feces	-	-
15	<i>Escherichia blattae</i>	ATCC 29907	Insect	-	-
16	<i>Escherichia hermannii</i>	ATCC 33650	Mouse Brain	-	-
17	<i>Escherichia vulneris</i>	ATCC 29943	Clinical Isolate, Human Wound	-	-
18	<i>Haemophilus influenzae</i>	ATCC 19418	Not Available	-	-
19	<i>Hafnia alvei</i>	ATCC 51815	Milk	-	-
20	<i>Klebsiella pneumoniae</i>	ATCC 4352	Not Available	-	-
21	<i>Kluyvera intermedia</i>	ATCC 33110	Surface Water	-	-
22	<i>Mycobacterium smegmatis</i>	ATCC 19420	Not Available	-	-
23	<i>Pantoea agglomerans</i>	ATCC 19552	Sewage	-	-
24	<i>Proteus mirabilis</i>	ATCC 7002	Urine	-	-
25	<i>Providencia rettgeri</i>	ATCC 14505	Not Available	-	-
26	<i>Serratia marcescens</i>	ATCC 13880	Not Available	-	-
27	<i>Shigella boydii</i>	ATCC 9290	Pork Liver	-	-
28	<i>Siccibacter paucimobilis</i>	CCUG ⁶ 54945	Not Available	-	-
29	<i>Sphingomonas paucimobilis</i>	ATCC 29837	Hospital Respirator	-	-
30	<i>Vibrio vulnificus</i>	QL 02111-1A	Seafood Product	-	-

1. ATCC – American Type Culture Collection, 2. MSU – Michigan State University Culture Collection, 3. NCTC – National Collection of Type Cultures, 4. PSU - Pennsylvania State University Culture Collection, 5. QL – Q Laboratories Culture Collection, 6. CCUG - Culture Collection University of Gothenburg

Table 13. Fresh Raw Ground Beef Results: Analyzed on the Agilent AriaMX by the BACGene Evaluation Sheet¹, Candidate vs. Reference – POD Results (1)

STEC Strain	PREraser Treatment	Enrichment Time Point	MPN ^a / Test Portion	N ^b	Candidate ²			Reference			dPOD _c ^f	95% CI ^g
					x ^c	POD _c ^d	95% CI	x ^c	POD _R ^e	95% CI		
<i>E. coli</i> O157:H7 (NM) NCTC ³ 13125	Without PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.39 (0.17, 0.68)	20	8	0.40	0.22, 0.61	7	0.35	0.18, 0.57	0.05	-0.23, 0.32
			4.12 (1.58, 10.72)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.39 (0.17, 0.68)	20	8	0.40	0.22, 0.61	7	0.35	0.18, 0.57	0.05	-0.23, 0.32
			4.12 (1.58, 10.72)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.39 (0.17, 0.68)	20	8	0.40	0.22, 0.61	7	0.35	0.18, 0.57	0.05	-0.23, 0.32
			4.12 (1.58, 10.72)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.39 (0.17, 0.68)	20	8	0.40	0.22, 0.61	7	0.35	0.18, 0.57	0.05	-0.23, 0.32
			4.12 (1.58, 10.72)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹The BACGene Evaluation Sheet and FastFinder software produced the same results.²The BACGene STEC workflow and the BACGene *E. coli* O157:H7 workflow produced the same results.³NCTC – National Collection of Type Cultures^aMPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval^bN = Number of test portions^cx = Number of positive test portions^dPOD_c = Candidate method confirmed positive outcomes divided by the total number of trials^ePOD_R = Reference method confirmed positive outcomes divided by the total number of trials^fdPOD_c = Difference between the confirmed candidate method result and reference method confirmed result POD values^g95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level**Table 14. Fresh Raw Ground Beef Results: Analyzed on the Agilent AriaMX by the BACGene Evaluation Sheet¹, Presumptive vs. Confirmed – POD Results (1)**

STEC Strain	PREraser Treatment	Enrichment Time Point	MPN ^a / Test Portion	N ^b	Presumptive ²			Confirmed ³			dPOD _{CP} ^f	95% CI ^g
					x ^c	POD _{CP} ^d	95% CI	x ^c	POD _{CC} ^e	95% CI		
<i>E. coli</i> O157:H7 (NM) NCTC ⁴ 13125	Without PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.39 (0.17, 0.68)	20	8	0.40	0.22, 0.61	8	0.40	0.22, 0.61	0.00	N/A
			4.12 (1.58, 10.72)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.39 (0.17, 0.68)	20	8	0.40	0.22, 0.61	8	0.40	0.22, 0.61	0.00	N/A
			4.12 (1.58, 10.72)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.39 (0.17, 0.68)	20	8	0.40	0.22, 0.61	8	0.40	0.22, 0.61	0.00	N/A
			4.12 (1.58, 10.72)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.39 (0.17, 0.68)	20	8	0.40	0.22, 0.61	8	0.40	0.22, 0.61	0.00	N/A
			4.12 (1.58, 10.72)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹The BACGene Evaluation Sheet and FastFinder software produced the same results.²The BACGene STEC workflow and the BACGene *E. coli* O157:H7 workflow produced the same results.³Alternative and traditional confirmation procedures produced the same results.⁴NCTC – National Collection of Type Cultures^aMPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval^bN = Number of test portions^cx = Number of positive test portions^dPOD_{CP} = Candidate method presumptive positive outcomes divided by the total number of trials^ePOD_{CC} = Candidate method confirmed positive outcomes divided by the total number of trials^fdPOD_{CP} = Difference between the candidate method presumptive result and candidate method confirmed result POD values^g95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level

Table 15. Fresh Raw Ground Beef Results: Analyzed on the Agilent AriaMX by the BACGene Evaluation Sheet¹, Candidate vs. Reference – POD Results (1)

STEC Strain	PRE _{Eraser} Treatment	Enrichment Time Point	MPN ^a / Test Portion	N ^b	Candidate			Reference			dPOD _c ^f	95% CI ^g
					x ^c	POD _c ^d	95% CI	x ^c	POD _R ^e	95% CI		
<i>E. coli</i> O26 MSU ² TW00971	Without PRE _{Eraser}	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.49 (0.25, 0.84)	20	7	0.35	0.18, 0.57	9	0.45	0.26, 0.66	-0.10	-0.37, 0.19
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.49 (0.25, 0.84)	20	6	0.30	0.15, 0.52	9	0.45	0.26, 0.66	-0.15	-0.41, 0.14
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PRE _{Eraser}	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.49 (0.25, 0.84)	20	7	0.35	0.18, 0.57	9	0.45	0.26, 0.66	-0.10	-0.37, 0.19
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.49 (0.25, 0.84)	20	7	0.35	0.18, 0.57	9	0.45	0.26, 0.66	-0.10	-0.37, 0.19
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹The BACGene Evaluation Sheet and FastFinder software produced the same results.²MSU – Michigan State University Culture Collection^aMPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval^bN = Number of test portions^cx = Number of positive test portions^dPOD_c = Candidate method confirmed positive outcomes divided by the total number of trials^ePOD_R = Reference method confirmed positive outcomes divided by the total number of trials^fdPOD_c = Difference between the confirmed candidate method result and reference method confirmed result POD values^g95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level**Table 16. Fresh Raw Ground Beef Results: Analyzed on the Agilent AriaMX by the BACGene Evaluation Sheet¹, Presumptive vs. Confirmed – POD Results (1)**

STEC Strain	PRE _{Eraser} Treatment	Enrichment Time Point	MPN ^a / Test Portion	N ^b	Presumptive			Confirmed ²			dPOD _{CP} ^f	95% CI ^g
					x ^c	POD _{CP} ^d	95% CI	x ^c	POD _{CC} ^e	95% CI		
<i>E. coli</i> O26 MSU ³ TW00971	Without PRE _{Eraser}	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.49 (0.25, 0.84)	20	7	0.35	0.18, 0.57	7	0.35	0.18, 0.57	0.00	N/A
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.49 (0.25, 0.84)	20	6	0.30	0.15, 0.52	7	0.35	0.18, 0.57	-0.05	-0.21, 0.11
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PRE _{Eraser}	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.49 (0.25, 0.84)	20	7	0.35	0.18, 0.57	7	0.35	0.18, 0.57	0.00	N/A
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.49 (0.25, 0.84)	20	7	0.35	0.18, 0.57	7	0.35	0.18, 0.57	0.00	N/A
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹The BACGene Evaluation Sheet and FastFinder software produced the same results.²Alternative and traditional confirmation procedures produced the same results.³MSU – Michigan State University Culture Collection^aMPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval^bN = Number of test portions^cx = Number of positive test portions^dPOD_{CP} = Candidate method presumptive positive outcomes divided by the total number of trials^ePOD_{CC} = Candidate method confirmed positive outcomes divided by the total number of trials^fdPOD_{CP} = Difference between the candidate method presumptive result and candidate method confirmed result POD values^g95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level

Table 17. Fresh Raw Beef Trim Results: Analyzed by the BACGene Evaluation Sheet¹, Candidate vs. Reference – POD Results (1)

STEC Strain	Thermocycler	Enrichment Time Point	MPN ^a / Test Portion	N ^b	Candidate ²			Reference			dPOD _c ^f	95% CI ^g
					x ^c	POD _c ^d	95% CI	x ^c	POD _R ^e	95% CI		
<i>E. coli</i> O157:H7 ATCC ³ 35150	Agilent AriaMX™	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.68 (0.38, 1.14)	20	9	0.45	0.26, 0.66	11	0.55	0.34, 0.74	-0.10	-0.37, 0.19
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.68 (0.38, 1.14)	20	8	0.40	0.22, 0.61	11	0.55	0.34, 0.74	-0.15	-0.41, 0.14
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	Bio-Rad CFX96 Touch™ Deep Well	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.68 (0.38, 1.14)	20	9	0.45	0.26, 0.66	11	0.55	0.34, 0.74	-0.10	-0.37, 0.19
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.68 (0.38, 1.14)	20	9	0.45	0.26, 0.66	11	0.55	0.34, 0.74	-0.10	-0.37, 0.19
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹The BACGene Evaluation Sheet and FastFinder software produced the same results.²The BACGene STEC workflow and the BACGene *E. coli* O157:H7 workflow produced the same results.³ ATCC – American Type Culture Collection^aMPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval^bN = Number of test portions^cx = Number of positive test portions^dPOD_c = Candidate method confirmed positive outcomes divided by the total number of trials^ePOD_R = Reference method confirmed positive outcomes divided by the total number of trials^fdPOD_c = Difference between the confirmed candidate method result and reference method confirmed result POD values^g95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level**Table 18. Fresh Raw Beef Trim Results: Analyzed by the BACGene Evaluation Sheet¹, Presumptive vs. Confirmed – POD Results (1)**

STEC Strain	Thermocycler	Enrichment Time Point	MPN ^a / Test Portion	N ^b	Presumptive ²			Confirmed ³			dPOD _{cp} ^f	95% CI ^g
					x ^c	POD _{cp} ^d	95% CI	x ^c	POD _{cc} ^e	95% CI		
<i>E. coli</i> O157:H7 ATCC ⁴ 35150	Agilent AriaMX™	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.68 (0.38, 1.14)	20	9	0.45	0.26, 0.66	9	0.45	0.26, 0.66	0.00	N/A
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.68 (0.38, 1.14)	20	8	0.40	0.22, 0.61	9	0.45	0.26, 0.66	-0.05	-0.21, 0.11
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	Bio-Rad CFX96 Touch™ Deep Well	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.68 (0.38, 1.14)	20	9	0.45	0.26, 0.66	9	0.45	0.26, 0.66	0.00	N/A
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.68 (0.38, 1.14)	20	9	0.45	0.26, 0.66	9	0.45	0.26, 0.66	0.00	N/A
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹The BACGene Evaluation Sheet and FastFinder software produced the same results.²The BACGene STEC workflow and the BACGene *E. coli* O157:H7 workflow produced the same results.³Alternative and traditional confirmation procedures produced the same results.⁴ ATCC – American Type Culture Collection^aMPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval^bN = Number of test portions^cx = Number of positive test portions^dPOD_{cp} = Candidate method presumptive positive outcomes divided by the total number of trials^ePOD_{cc} = Candidate method confirmed positive outcomes divided by the total number of trials^fdPOD_{cp} = Difference between the candidate method presumptive result and candidate method confirmed result POD values^g95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level

Table 19. Fresh Raw Beef Trim Results: Analyzed by the BACGene Evaluation Sheet¹, Candidate vs. Reference – POD Results (1)

STEC Strain	Thermocycler	Enrichment Time Point	MPN ^a / Test Portion	N ^b	Candidate			Reference			dPOD _c ^f	95% CI ^g
					x ^c	POD _c ^d	95% CI	x ^c	POD _R ^e	95% CI		
<i>E. coli</i> O45 MSU ² TW10121	Agilent AriaMX™	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.81 (0.47, 1.36)	20	7	0.35	0.18, 0.57	12	0.60	0.39, 0.78	-0.25	-0.50, 0.05
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.81 (0.47, 1.36)	20	7	0.35	0.18, 0.57	12	0.60	0.39, 0.78	-0.25	-0.50, 0.05
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	Bio-Rad CFX96 Touch™ Deep Well	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.81 (0.47, 1.36)	20	7	0.35	0.18, 0.57	12	0.60	0.39, 0.78	-0.25	-0.50, 0.05
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.81 (0.47, 1.36)	20	7	0.35	0.18, 0.57	12	0.60	0.39, 0.78	-0.25	-0.50, 0.05
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹The BACGene Evaluation Sheet and FastFinder software produced the same results.²MSU - Michigan State Culture Collection^aMPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval^bN = Number of test portions^cx = Number of positive test portions^dPOD_c = Candidate method confirmed positive outcomes divided by the total number of trials^ePOD_R = Reference method confirmed positive outcomes divided by the total number of trials^fdPOD_c = Difference between the confirmed candidate method result and reference method confirmed result POD values^g95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level**Table 20. Fresh Raw Beef Trim Results: Analyzed by the BACGene Evaluation Sheet¹, Presumptive vs. Confirmed – POD Results (1)**

STEC Strain	Thermocycler	Enrichment Time Point	MPN ^a / Test Portion	N ^b	Presumptive			Confirmed ²			dPOD _{CP} ^f	95% CI ^g
					x ^c	POD _{CP} ^d	95% CI	x ^c	POD _{CC} ^e	95% CI		
<i>E. coli</i> O45 MSU ³ TW10121	Agilent AriaMX™	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.81 (0.47, 1.36)	20	7	0.35	0.18, 0.57	7	0.35	0.18, 0.57	0.00	N/A
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.81 (0.47, 1.36)	20	7	0.35	0.18, 0.57	7	0.35	0.18, 0.57	0.00	N/A
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	Bio-Rad CFX96 Touch™ Deep Well	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.81 (0.47, 1.36)	20	7	0.35	0.18, 0.57	7	0.35	0.18, 0.57	0.00	N/A
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.81 (0.47, 1.36)	20	7	0.35	0.18, 0.57	7	0.35	0.18, 0.57	0.00	N/A
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹The BACGene Evaluation Sheet and FastFinder software produced the same results.²Alternative and traditional confirmation procedures produced the same results.³MSU - Michigan State Culture Collection^aMPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval^bN = Number of test portions^cx = Number of positive test portions^dPOD_{CP} = Candidate method presumptive positive outcomes divided by the total number of trials^ePOD_{CC} = Candidate method confirmed positive outcomes divided by the total number of trials^fdPOD_{CP} = Difference between the candidate method presumptive result and candidate method confirmed result POD values^g95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level

Table 21. Beef Carcass Sampling Cloth Results: Analyzed on the Agilent AriaMX by the BACGene Evaluation Sheet¹, Candidate vs. Reference – POD Results (1)

Matrix	STEC Strain	Enrichment Time Point	CFU/Test Area	N ^b	Candidate			Reference			dPOD _c ^f	95% CI ^g
					x ^c	POD _c ^d	95% CI	x ^c	POD _R ^e	95% CI		
Beef Carcass Sampling Cloths	<i>E. coli</i> O157:H7 USDA 465-97	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			56	20	9	0.45	0.26, 0.66	7	0.35	0.18, 0.57	0.10	-0.19, 0.37
			480	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			56	20	9	0.45	0.26, 0.66	7	0.35	0.18, 0.57	0.10	-0.19, 0.37
			480	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	<i>E. coli</i> O111 MSU ² DEC 8D ³	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			51	20	10	0.50	0.30, 0.70	6	0.30	0.15, 0.52	0.20	-0.10, 0.45
			530	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			51	20	10	0.50	0.30, 0.70	6	0.30	0.15, 0.52	0.20	-0.10, 0.45
			530	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹The BACGene Evaluation Sheet and FastFinder software produced the same results.²MSU – Michigan State Culture Collection³The BACGene STEC workflow and the BACGene *E. coli* O157:H7 workflow produced the same results.⁴MPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval⁵N = Number of test portions⁶x = Number of positive test portions⁷POD_c = Candidate method confirmed positive outcomes divided by the total number of trials⁸POD_R = Reference method confirmed positive outcomes divided by the total number of trials⁹dPOD_c = Difference between the confirmed candidate method result and reference method confirmed result POD values¹⁰95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level**Table 22. Beef Carcass Sampling Cloth Results: Analyzed on the Agilent AriaMX by the BACGene Evaluation Sheet¹, Presumptive vs. Confirmed – POD Results (1)**

Matrix	STEC Strain	Enrichment Time Point	CFU/Test Area	N ^b	Presumptive			Confirmed ²			dPOD _{cp} ^f	95% CI ^g
					x ^c	POD _{cp} ^d	95% CI	x ^c	POD _{cc} ^e	95% CI		
Beef Carcass Sampling Cloths	<i>E. coli</i> O157:H7 USDA 465-97	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			56	20	9	0.45	0.26, 0.66	9	0.45	0.26, 0.66	0.00	N/A
			480	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			56	20	9	0.45	0.26, 0.66	9	0.45	0.26, 0.66	0.00	N/A
			480	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	<i>E. coli</i> O111 MSU ³ DEC 8D ⁴	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			51	20	10	0.50	0.30, 0.70	10	0.50	0.30, 0.70	0.00	N/A
			530	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			51	20	10	0.50	0.30, 0.70	10	0.50	0.30, 0.70	0.00	N/A
			530	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹The BACGene Evaluation Sheet and FastFinder software produced the same results.² Alternative and traditional confirmation procedures produced the same results.³MSU – Michigan State Culture Collection⁴The BACGene STEC workflow and the BACGene *E. coli* O157:H7 workflow produced the same results.⁵MPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval⁶N = Number of test portions⁷x = Number of positive test portions⁸POD_{cp} = Candidate method presumptive positive outcomes divided by the total number of trials⁹POD_{cc} = Candidate method confirmed positive outcomes divided by the total number of trials¹⁰dPOD_{cp} = Difference between the candidate method presumptive result and candidate method confirmed result POD values¹¹95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level

Table 23. Romaine Lettuce Results: Analyzed by the BACGene Evaluation Sheet¹, Candidate vs. Reference – POD Results (1)

STEC Strain	Thermocycler	Enrichment Time Point	MPN ^a / Test Portion	N ^b	Candidate ²			Reference			dPOD ^c	95% CI ^g
					x ^c	POD ^d	95% CI	x ^c	POD ^e	95% CI		
<i>E. coli</i> O157:H7 ATCC ³ 43895	Agilent AriaMX™	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.45 (0.22, 0.78)	20	7	0.35	0.18, 0.57	7	0.35	0.18, 0.57	0.00	N/A
			3.70 (1.52, 9.02)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.45 (0.22, 0.78)	20	7	0.35	0.18, 0.57	7	0.35	0.18, 0.57	0.00	N/A
			3.70 (1.52, 9.02)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	Bio-Rad CFX96 Touch™ Deep Well	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.45 (0.22, 0.78)	20	7	0.35	0.18, 0.57	7	0.35	0.18, 0.57	0.00	N/A
			3.70 (1.52, 9.02)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.45 (0.22, 0.78)	20	7	0.35	0.18, 0.57	7	0.35	0.18, 0.57	0.00	N/A
			3.70 (1.52, 9.02)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹The BACGene Evaluation Sheet and FastFinder software produced the same results.²The BACGene STEC workflow and the BACGene *E. coli* O157:H7 workflow produced the same results.³ATCC – American Type Culture Collection^aMPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval^bN = Number of test portions^cx = Number of positive test portions^dPOD_c = Candidate method confirmed positive outcomes divided by the total number of trials^ePOD_R = Reference method confirmed positive outcomes divided by the total number of trials^fdPOD_c = Difference between the confirmed candidate method result and reference method confirmed result POD values^g95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level**Table 24. Romaine Lettuce Results: Analyzed by the BACGene Evaluation Sheet¹, Presumptive vs. Confirmed – POD Results (1)**

STEC Strain	Thermocycler	Enrichment Time Point	MPN ^a / Test Portion	N ^b	Presumptive ²			Confirmed ³			dPOD _{CP} ^f	95% CI ^g
					x ^c	POD _{CP} ^d	95% CI	x ^c	POD _{CC} ^e	95% CI		
<i>E. coli</i> O157:H7 ATCC ⁴ 43895	Agilent AriaMX™	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.45 (0.22, 0.78)	20	7	0.35	0.18, 0.57	7	0.35	0.18, 0.57	0.00	N/A
			3.70 (1.52, 9.02)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.45 (0.22, 0.78)	20	7	0.35	0.18, 0.57	7	0.35	0.18, 0.57	0.00	N/A
			3.70 (1.52, 9.02)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	Bio-Rad CFX96 Touch™ Deep Well	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.45 (0.22, 0.78)	20	7	0.35	0.18, 0.57	7	0.35	0.18, 0.57	0.00	N/A
			3.70 (1.52, 9.02)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.45 (0.22, 0.78)	20	7	0.35	0.18, 0.57	7	0.35	0.18, 0.57	0.00	N/A
			3.70 (1.52, 9.02)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹The BACGene Evaluation Sheet and FastFinder software produced the same results.²The BACGene STEC workflow and the BACGene *E. coli* O157:H7 workflow produced the same results.³Alternative and traditional confirmation procedures produced the same results.⁴ATCC – American Type Culture Collection^aMPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval^bN = Number of test portions^cx = Number of positive test portions^dPOD_{CP} = Candidate method presumptive positive outcomes divided by the total number of trials^ePOD_{CC} = Candidate method confirmed positive outcomes divided by the total number of trials^fdPOD_{CP} = Difference between the candidate method presumptive result and candidate method confirmed result POD values^g95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level

Table 25. Romaine Lettuce Results: Analyzed by the BACGene Evaluation Sheet, Candidate vs. Reference – POD Results (1)

STEC Strain	Thermocycler	Enrichment Time Point	MPN ^a / Test Portion	N ^b	Candidate			Reference			dPOD ^c	95% CI ^e
					x ^c	POD _C ^d	95% CI	x ^c	POD _R ^e	95% CI		
<i>E. coli</i> O121 PSU ¹ 5.0959	Agilent AriaMX™	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.54 (0.29, 0.92)	20	9	0.45	0.26, 0.66	7	0.35	0.18, 0.57	0.10	-0.19, 0.37
			3.70 (1.52, 9.02)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.54 (0.29, 0.92)	20	10	0.50	0.30, 0.70	7	0.35	0.18, 0.57	0.15	-0.15, 0.41
			3.70 (1.52, 9.02)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	Bio-Rad CFX96 Touch™ Deep Well	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.54 (0.29, 0.92)	20	10	0.50	0.30, 0.70	7	0.35	0.18, 0.57	0.15	-0.15, 0.41
			3.70 (1.52, 9.02)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.54 (0.29, 0.92)	20	10	0.50	0.30, 0.70	7	0.35	0.18, 0.57	0.15	-0.15, 0.41
			3.70 (1.52, 9.02)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹PSU – Pennsylvania State university Culture Collection^aMPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval^bN = Number of test portions^cx = Number of positive test portions^dPOD_C = Candidate method confirmed positive outcomes divided by the total number of trials^ePOD_R = Reference method confirmed positive outcomes divided by the total number of trials^fdPOD_C = Difference between the confirmed candidate method result and reference method confirmed result POD values^g95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level**Table 26. Romaine Lettuce Results: Analyzed by the FastFinder Software, Candidate vs. Reference – POD Results (1)**

STEC Strain	Thermocycler	Enrichment Time Point	MPN ^a / Test Portion	N ^b	Candidate			Reference			dPOD ^c	95% CI ^e
					x ^c	POD _C ^d	95% CI	x ^c	POD _R ^e	95% CI		
<i>E. coli</i> O121 PSU ¹ 5.0959	Agilent AriaMX™	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.54 (0.29, 0.92)	20	10	0.50	0.30, 0.70	7	0.35	0.18, 0.57	0.15	-0.15, 0.41
			3.70 (1.52, 9.02)	5	4	8.00	0.38, 1.00	5	1.00	0.57, 1.00	-0.20	-0.62, 0.28
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.54 (0.29, 0.92)	20	10	0.50	0.30, 0.70	7	0.35	0.18, 0.57	0.15	-0.15, 0.41
			3.70 (1.52, 9.02)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	Bio-Rad CFX96 Touch™ Deep Well	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.54 (0.29, 0.92)	20	10	0.50	0.30, 0.70	7	0.35	0.18, 0.57	0.15	-0.15, 0.41
			3.70 (1.52, 9.02)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.54 (0.29, 0.92)	20	10	0.50	0.30, 0.70	7	0.35	0.18, 0.57	0.15	-0.15, 0.41
			3.70 (1.52, 9.02)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹PSU – Pennsylvania State university Culture Collection^aMPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval^bN = Number of test portions^cx = Number of positive test portions^dPOD_C = Candidate method confirmed positive outcomes divided by the total number of trials^ePOD_R = Reference method confirmed positive outcomes divided by the total number of trials^fdPOD_C = Difference between the confirmed candidate method result and reference method confirmed result POD values^g95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level

Table 27. Romaine Lettuce Results: Analyzed by the BACGene Evaluation Sheet, Presumptive vs. Confirmed – POD Results (1)

STEC Strain	Thermocycler	Enrichment Time Point	MPN ^a / Test Portion	N ^b	Presumptive			Confirmed ¹			dPOD _{CP} ^f	95% CI ^g
					x ^c	POD _{CP} ^d	95% CI	x ^c	POD _{CC} ^e	95% CI		
<i>E. coli</i> O121 PSU ² 5.0959	Agilent AriaMX™	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.54 (0.29, 0.92)	20	12	0.60	0.39, 0.78	10	0.50	0.30, 0.70	0.10	-0.13, 0.33
			3.70 (1.52, 9.02)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.54 (0.29, 0.92)	20	12	0.60	0.39, 0.78	10	0.50	0.30, 0.70	0.10	-0.08, 0.28
			3.70 (1.52, 9.02)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	Bio-Rad CFX96 Touch™ Deep Well	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.54 (0.29, 0.92)	20	11	0.55	0.34, 0.74	10	0.50	0.30, 0.70	0.05	-0.11, 0.21
			3.70 (1.52, 9.02)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.54 (0.29, 0.92)	20	10	0.50	0.30, 0.70	10	0.50	0.30, 0.70	0.00	N/A
			3.70 (1.52, 9.02)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹Alternative and traditional confirmation procedures produced the same results.²PSU – Pennsylvania State university Culture Collection^aMPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval^bN = Number of test portions^cx = Number of positive test portions^dPOD_{CP} = Candidate method presumptive positive outcomes divided by the total number of trials^ePOD_{CC} = Candidate method confirmed positive outcomes divided by the total number of trials^fdPOD_{CP} = Difference between the candidate method presumptive result and candidate method confirmed result POD values^g95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level**Table 28. Romaine Lettuce Results: Analyzed by the FastFinder Software, Presumptive vs. Confirmed – POD Results (1)**

STEC Strain	Thermocycler	Enrichment Time Point	MPN ^a / Test Portion	N ^b	Presumptive			Confirmed ¹			dPOD _{CP} ^f	95% CI ^g
					x ^c	POD _{CP} ^d	95% CI	x ^c	POD _{CC} ^e	95% CI		
<i>E. coli</i> O121 PSU ² 5.0959	Agilent AriaMX™	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.54 (0.29, 0.92)	20	11	0.55	0.34, 0.77	10	0.50	0.30, 0.70	0.05	-0.11, 0.21
			3.70 (1.52, 9.02)	5	4	8.00	0.38, 1.00	5	1.00	0.57, 1.00	-0.20	-0.76, 0.36
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.54 (0.29, 0.92)	20	10	0.50	0.30, 0.70	10	0.50	0.30, 0.70	0.00	N/A
			3.70 (1.52, 9.02)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	Bio-Rad CFX96 Touch™ Deep Well	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.54 (0.29, 0.92)	20	11	0.55	0.34, 0.74	10	0.50	0.30, 0.70	0.05	-0.11, 0.21
			3.70 (1.52, 9.02)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.54 (0.29, 0.92)	20	10	0.50	0.30, 0.70	10	0.50	0.30, 0.70	0.00	N/A
			3.70 (1.52, 9.02)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹Alternative and traditional confirmation procedures produced the same results.²PSU – Pennsylvania State university Culture Collection^aMPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval^bN = Number of test portions^cx = Number of positive test portions^dPOD_{CP} = Candidate method presumptive positive outcomes divided by the total number of trials^ePOD_{CC} = Candidate method confirmed positive outcomes divided by the total number of trials^fdPOD_{CP} = Difference between the candidate method presumptive result and candidate method confirmed result POD values^g95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level

Table 29. All Purpose Flour Results: Analyzed on the Bio-Rad CFX96 Touch™ Deep Well by the BACGene Evaluation Sheet¹, Candidate vs. Reference – POD Results (1)

STEC Strain	PREraser Treatment	Enrichment Time Point	MPN ^a / Test Portion	N ^b	Candidate ²			Reference			dPOD _c ^f	95% CI ^g
					x ^c	POD _c ^d	95% CI	x ^c	POD _R ^e	95% CI		
<i>E. coli</i> O157:H7 ATCC ³ BAA-460	Without PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.34 (0.14, 0.61)	20	10	0.50	0.30, 0.70	6	0.30	0.15, 0.52	0.20	-0.10, 0.45
			4.12 (1.58, 10.72)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.34 (0.14, 0.61)	20	10	0.50	0.30, 0.70	6	0.30	0.15, 0.52	0.20	-0.10, 0.45
			4.12 (1.58, 10.72)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	-0.43, 0.43
			0.34 (0.14, 0.61)	20	10	0.50	0.30, 0.70	6	0.30	0.15, 0.52	0.20	-0.10, 0.45
			4.12 (1.58, 10.72)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.34 (0.14, 0.61)	20	10	0.50	0.30, 0.70	6	0.30	0.15, 0.52	0.20	-0.10, 0.45
			4.12 (1.58, 10.72)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹The BACGene Evaluation Sheet and FastFinder software produced the same results.²The BACGene STEC workflow and the BACGene *E. coli* O157:H7 workflow produced the same results.³ATCC – American Type Culture Collection^aMPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval^bN = Number of test portions^cx = Number of positive test portions^dPOD_c = Candidate method confirmed positive outcomes divided by the total number of trials^ePOD_R = Reference method confirmed positive outcomes divided by the total number of trials^fdPOD_c = Difference between the confirmed candidate method result and reference method confirmed result POD values^g95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level**Table 30. All Purpose Flour Results: Analyzed on the Bio-Rad CFX96 Touch™ Deep Well by the BACGene Evaluation Sheet, Presumptive vs. Confirmed – POD Results (1)**

STEC Strain	PREraser Treatment	Enrichment Time Point	MPN ^a / Test Portion	N ^b	Presumptive ¹			Confirmed ²			dPOD _{CP} ^f	95% CI ^g
					x ^c	POD _{CP} ^d	95% CI	x ^c	POD _{CC} ^e	95% CI		
<i>E. coli</i> O157:H7 ATCC ³ BAA-460	Without PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.34 (0.14, 0.61)	20	10	0.50	0.30, 0.70	10	0.50	0.30, 0.70	0.00	N/A
			4.12 (1.58, 10.72)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.34 (0.14, 0.61)	20	10	0.50	0.30, 0.77	10	0.50	0.30, 0.70	0.00	N/A
			4.12 (1.58, 10.72)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.34 (0.14, 0.61)	20	10	0.50	0.30, 0.70	10	0.50	0.30, 0.70	0.00	N/A

			4.12 (1.58, 10.72)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.34 (0.14, 0.61)	20	11	0.55	0.34, 0.74	10	0.50	0.30, 0.70	0.05	N/A
			4.12 (1.58, 10.72)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹The BACGene STEC workflow and the BACGene *E. coli* O157:H7 workflow produced the same results.

²Alternative and traditional confirmation procedures produced the same results.

³ATCC – American Type Culture Collection

^aMPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval

^bN = Number of test portions

^cx = Number of positive test portions

^dPOD_{CP} = Candidate method presumptive positive outcomes divided by the total number of trials

^ePOD_{CC} = Candidate method confirmed positive outcomes divided by the total number of trials

^fdPOD_{CP} = Difference between the candidate method presumptive result and candidate method confirmed result POD values

^g95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level

Table 31. All Purpose Flour Results: Analyzed on the Bio-Rad CFX96 Touch Deep Well by the FastFinder Software, Presumptive vs. Confirmed – POD Results (1)

STEC Strain	PREraser Treatment	Enrichment Time Point	MPN ^a / Test Portion	N ^b	Presumptive ¹			Confirmed ²			dPOD _{CP} ^f	95% CI ^g
					x ^c	POD _{CP} ^d	95% CI	x ^c	POD _{CC} ^e	95% CI		
<i>E. coli</i> O157:H7 ATCC ³ BAA-460	Without PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.34 (0.14, 0.61)	20	10	0.50	0.30, 0.70	10	0.50	0.30, 0.70	0.00	N/A
			4.12 (1.58, 10.72)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.34 (0.14, 0.61)	20	10	0.50	0.30, 0.70	10	0.50	0.30, 0.70	0.00	N/A
			4.12 (1.58, 10.72)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.34 (0.14, 0.61)	20	10	0.50	0.30, 0.70	10	0.50	0.30, 0.70	0.00	N/A
			4.12 (1.58, 10.72)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.34 (0.14, 0.61)	20	11	0.55	0.34, 0.74	10	0.50	0.30, 0.70	0.05	N/A
			4.12 (1.58, 10.72)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹The BACGene STEC workflow and the BACGene *E. coli* O157:H7 workflow produced the same results.

²Alternative and traditional confirmation procedures produced the same results.

³ATCC – American Type Culture Collection

^aMPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval

^bN = Number of test portions

^cx = Number of positive test portions

^dPOD_{CP} = Candidate method presumptive positive outcomes divided by the total number of trials

^ePOD_{CC} = Candidate method confirmed positive outcomes divided by the total number of trials

^fdPOD_{CP} = Difference between the candidate method presumptive result and candidate method confirmed result POD values

^g95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level

Table 32. All Purpose Flour Results: Analyzed on the Bio-Rad CFX96 Touch Deep Well by the BACGene Evaluation Sheet¹, Candidate vs. Reference – POD Results (1)

STEC Strain	PREraser Treatment	Enrichment Time Point	MPN ^a / Test Portion	N ^b	Candidate			Reference			dPOD _c ^f	95% CI ^g
					x ^c	POD _c ^d	95% CI	x ^c	POD _R ^e	95% CI		
<i>E. coli</i> O103 MSU ² TW09101	Without PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.40 (0.18, 0.69)	20	7	0.35	0.18, 0.57	6	0.30	0.15, 0.52	0.05	-0.23, 0.32
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.40 (0.18, 0.69)	20	7	0.35	0.18, 0.57	6	0.30	0.15, 0.52	0.05	-0.23, 0.32
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.40 (0.18, 0.69)	20	7	0.35	0.18, 0.57	6	0.30	0.15, 0.52	0.05	-0.23, 0.32
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.40 (0.18, 0.69)	20	7	0.35	0.18, 0.57	6	0.30	0.15, 0.52	0.05	-0.23, 0.32
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹The BACGene Evaluation Sheet and FastFinder software produced the same results.²MSU – Michigan State Culture Collection^aMPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval^bN = Number of test portions^cx = Number of positive test portions^dPOD_c = Candidate method confirmed positive outcomes divided by the total number of trials^ePOD_R = Reference method confirmed positive outcomes divided by the total number of trials^fdPOD_c = Difference between the confirmed candidate method result and reference method confirmed result POD values^g95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level**Table 33. All Purpose Flour Results: Analyzed on the Bio-Rad CFX96 Touch Deep Well by the BACGene Evaluation Sheet¹, Presumptive vs. Confirmed – POD Results (1)**

STEC Strain	PREraser Treatment	Enrichment Time Point	MPN ^a / Test Portion	N ^b	Presumptive			Confirmed ²			dPOD _{cp} ^f	95% CI ^g
					x ^c	POD _{cp} ^d	95% CI	x ^c	POD _R ^e	95% CI		
<i>E. coli</i> O103 MSU ³ TW09101	Without PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.40 (0.18, 0.69)	20	8	0.40	0.22, 0.61	7	0.35	0.18, 0.57	0.05	-0.11, 0.21
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.40 (0.18, 0.69)	20	7	0.35	0.18, 0.57	7	0.35	0.18, 0.57	0.00	N/A
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PREraser ³	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.40 (0.18, 0.69)	20	7	0.35	0.18, 0.57	7	0.35	0.18, 0.57	0.00	N/A
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			0.40 (0.18, 0.69)	20	7	0.35	0.18, 0.57	7	0.35	0.18, 0.57	0.00	N/A
			6.53 (2.30, 18.54)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹The BACGene Evaluation Sheet and FastFinder software produced the same results.

²Alternative and traditional confirmation procedures produced the same results.

³MSU – Michigan State Culture Collection

^aMPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval

^bN = Number of test portions

^cx = Number of positive test portions

^dPOD_c = Candidate method confirmed positive outcomes divided by the total number of trials

^ePOD_r = Reference method confirmed positive outcomes divided by the total number of trials

^fdPOD_c = Difference between the confirmed candidate method result and reference method confirmed result POD values

^g95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level

Table 34. Stainless Steel Environmental Sponge Results: Analyzed on the Bio-Rad CFX96 Touch Deep Well by the BACGene Evaluation Sheet¹, Candidate vs. Reference – POD Results (1)

STEC Strain	PREraser Treatment	Enrichment Time Point	CFU/Test Area	N ^b	Candidate ²			Reference			dPOD _c ^f	95% CI ^g
					x ^c	POD _c ^d	95% CI	x ^c	POD _r ^e	95% CI		
<i>E. coli</i> O157:H7 ATCC ³ 43895 with 10x <i>C. freundii</i> ATCC 8090	Without PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			55 (<i>E. coli</i>) & 790 (<i>C. freundii</i>)	20	8	0.40	0.22, 0.61	8	0.40	0.22, 0.61	0.00	N/A
			400 (<i>E. coli</i>) & 7300 (<i>C. freundii</i>)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			55 (<i>E. coli</i>) & 790 (<i>C. freundii</i>)	20	8	0.40	0.22, 0.61	8	0.40	0.22, 0.61	0.00	N/A
			400 (<i>E. coli</i>) & 7300 (<i>C. freundii</i>)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			55 (<i>E. coli</i>) & 790 (<i>C. freundii</i>)	20	8	0.40	0.22, 0.61	8	0.40	0.22, 0.61	0.00	N/A
			400 (<i>E. coli</i>) & 7300 (<i>C. freundii</i>)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			55 (<i>E. coli</i>) & 790 (<i>C. freundii</i>)	20	8	0.40	0.22, 0.61	8	0.40	0.22, 0.61	0.00	N/A
			400 (<i>E. coli</i>) & 7300 (<i>C. freundii</i>)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹The BACGene Evaluation Sheet and FastFinder software produced the same results.

²The BACGene STEC workflow and the BACGene *E. coli* O157:H7 workflow produced the same results.

³ATCC – American Type Culture Collection

^aMPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval

^bN = Number of test portions

^cx = Number of positive test portions

^dPOD_c = Candidate method confirmed positive outcomes divided by the total number of trials

^ePOD_r = Reference method confirmed positive outcomes divided by the total number of trials

^fdPOD_c = Difference between the confirmed candidate method result and reference method confirmed result POD values

^g95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level

Table 35. Stainless Steel Environmental Sponge Results: Analyzed on the Bio-Rad CFX96 Touch Deep Well by the BACGene Evaluation Sheet¹, Presumptive vs. Confirmed – POD Results (1)

STEC Strain	PREraser Treatment	Enrichment Time Point	CFU/Test Area	N ^b	Presumptive ²			Confirmed ³			dPOD _{CP} ^f	95% CI ^g
					x ^c	POD _{CP} ^d	95% CI	x ^c	POD _{CC} ^e	95% CI		
<i>E. coli</i> O157:H7 ATCC ⁴ 43895 with 10x <i>C. freundii</i> ATCC 8090	Without PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			55 (<i>E. coli</i>) & 790 (<i>C. freundii</i>)	20	8	0.40	0.22, 0.61	8	0.40	0.22, 0.61	0.00	N/A
			400 (<i>E. coli</i>) & 7300 (<i>C. freundii</i>)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			55 (<i>E. coli</i>) & 790 (<i>C. freundii</i>)	20	8	0.40	0.22, 0.61	8	0.40	0.22, 0.61	0.00	N/A
			400 (<i>E. coli</i>) & 7300 (<i>C. freundii</i>)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			55 (<i>E. coli</i>) & 790 (<i>C. freundii</i>)	20	8	0.40	0.22, 0.61	8	0.40	0.22, 0.61	0.00	N/A
			400 (<i>E. coli</i>) & 7300 (<i>C. freundii</i>)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			55 (<i>E. coli</i>) & 790 (<i>C. freundii</i>)	20	8	0.40	0.22, 0.61	8	0.40	0.22, 0.61	0.00	N/A
			400 (<i>E. coli</i>) & 7300 (<i>C. freundii</i>)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹The BACGene Evaluation Sheet and FastFinder software produced the same results.²The BACGene STEC workflow and the BACGene *E. coli* O157:H7 workflow produced the same results.³Alternative and traditional confirmation procedures produced the same results.⁴ATCC – American Type Culture Collection^aMPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval^bN = Number of test portions^cx = Number of positive test portions^dPOD_{CP} = Candidate method presumptive positive outcomes divided by the total number of trials^ePOD_{CC} = Candidate method confirmed positive outcomes divided by the total number of trials^fdPOD_{CP} = Difference between the candidate method presumptive result and candidate method confirmed result POD values^g95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level**Table 36. Stainless Steel Environmental Sponge Results: Analyzed on the Bio-Rad CFX96 Touch Deep Well by the BACGene Evaluation Sheet¹, Candidate vs. Reference – POD Results (1)**

STEC Strain	PREraser Treatment	Enrichment Time Point	CFU/Test Area	N ^b	Candidate			Reference			dPOD _C ^f	95% CI ^g
					x ^c	POD _C ^d	95% CI	x ^c	POD _R ^e	95% CI		
<i>E. coli</i> O111 MSU ² TW06315 with 10x <i>C. freundii</i> ATCC 8090	Without PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			62 (<i>E. coli</i>) & 790 (<i>C. freundii</i>)	20	9	0.45	0.26, 0.66	7	0.35	0.18, 0.57	0.10	-0.19, 0.37
			440 (<i>E. coli</i>) & 7300 (<i>C. freundii</i>)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A

			62 (<i>E. coli</i>) & 790 (<i>C. freundii</i>)	20	9	0.45	0.26, 0.66	7	0.35	0.18, 0.57	0.10	-0.19, 0.37
			440 (<i>E. coli</i>) & 7300 (<i>C. freundii</i>)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			62 (<i>E. coli</i>) & 790 (<i>C. freundii</i>)	20	9	0.45	0.26, 0.66	7	0.35	0.18, 0.57	0.10	-0.19, 0.37
			440 (<i>E. coli</i>) & 7300 (<i>C. freundii</i>)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			62 (<i>E. coli</i>) & 790 (<i>C. freundii</i>)	20	9	0.45	0.26, 0.66	7	0.35	0.18, 0.57	0.10	-0.19, 0.37
			440 (<i>E. coli</i>) & 7300 (<i>C. freundii</i>)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹The BACGene Evaluation Sheet and FastFinder software produced the same results.²MSU – Michigan State Culture Collection³MPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval⁴N = Number of test portions⁵x = Number of positive test portions⁶POD_c = Candidate method confirmed positive outcomes divided by the total number of trials⁷POD_r = Reference method confirmed positive outcomes divided by the total number of trials⁸dPOD_c = Difference between the confirmed candidate method result and reference method confirmed result POD values⁹95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level**Table 37. Stainless Steel Environmental Sponge Results: Analyzed on the Bio-Rad CFX96 Touch Deep Well by the BACGene Evaluation Sheet, Presumptive vs. Confirmed – POD Results (1)**

STEC Strain	PREraser Treatment	Enrichment Time Point	CFU/ Test Area	N ^b	Presumptive			Confirmed ¹			dPOD _{CP} ^f	95% CI ^g
					x ^c	POD _{CP} ^d	95% CI	x ^c	POD _{CC} ^e	95% CI		
<i>E. coli</i> O111 MSU ² TW06315 with 10x <i>C. freundii</i> ATCC 8090	Without PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			62 (<i>E. coli</i>) & 790 (<i>C. freundii</i>)	20	9	0.45	0.26, 0.66	9	0.45	0.26, 0.66	0.00	N/A
			440 (<i>E. coli</i>) & 7300 (<i>C. freundii</i>)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			62 (<i>E. coli</i>) & 790 (<i>C. freundii</i>)	20	9	0.45	0.26, 0.66	9	0.45	0.26, 0.66	0.00	N/A
			440 (<i>E. coli</i>) & 7300 (<i>C. freundii</i>)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			62 (<i>E. coli</i>) & 790 (<i>C. freundii</i>)	20	9	0.45	0.26, 0.66	9	0.45	0.26, 0.66	0.00	N/A
			440 (<i>E. coli</i>) & 7300 (<i>C. freundii</i>)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			62 (<i>E. coli</i>) & 790 (<i>C. freundii</i>)	20	10	0.50	0.30, 0.70	9	0.45	0.26, 0.66	0.05	-0.11, 0.21
			440 (<i>E. coli</i>) & 7300 (<i>C. freundii</i>)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹Alternative and traditional confirmation procedures produced the same results.

²MSU – Michigan State Culture Collection

³MPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval

⁴N = Number of test portions

⁵x = Number of positive test portions

⁶POD_{CP} = Candidate method presumptive positive outcomes divided by the total number of trials

⁷POD_{CC} = Candidate method confirmed positive outcomes divided by the total number of trials

⁸dPOD_{CP} = Difference between the candidate method presumptive result and candidate method confirmed result POD values

⁹95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level

Table 38. Stainless Steel Environmental Sponge Results: Analyzed on the Bio-Rad CFX96 Touch Deep Well by the FastFinder Software, Presumptive vs. Confirmed – POD Results (1)

STEC Strain	PREraser Treatment	Enrichment Time Point	CFU/ Test Area	N ^b	Presumptive			Confirmed ¹			dPOD _{CP} ^f	95% CI ^g
					x ^c	POD _{CP} ^d	95% CI	x ^c	POD _{CC} ^e	95% CI		
<i>E. coli</i> O111 MSU ² TW06315 with 10x <i>C. freundii</i> ATCC 8090	Without PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			62 (<i>E. coli</i>) & 790 (<i>C. freundii</i>)	20	9	0.45	0.26, 0.66	9	0.45	0.26, 0.66	0.00	N/A
			440 (<i>E. coli</i>) & 7300 (<i>C. freundii</i>)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			62 (<i>E. coli</i>) & 790 (<i>C. freundii</i>)	20	9	0.45	0.26, 0.66	9	0.45	0.26, 0.66	0.00	N/A
			440 (<i>E. coli</i>) & 7300 (<i>C. freundii</i>)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			62 (<i>E. coli</i>) & 790 (<i>C. freundii</i>)	20	9	0.45	0.26, 0.66	9	0.45	0.26, 0.66	0.00	N/A
			440 (<i>E. coli</i>) & 7300 (<i>C. freundii</i>)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			62 (<i>E. coli</i>) & 790 (<i>C. freundii</i>)	20	11	0.55	0.34, 0.74	9	0.45	0.26, 0.66	0.10	-0.08, 0.28
			440 (<i>E. coli</i>) & 7300 (<i>C. freundii</i>)	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹Alternative and traditional confirmation procedures produced the same results.

²MSU – Michigan State Culture Collection

³MPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval

⁴N = Number of test portions

⁵x = Number of positive test portions

⁶POD_{CP} = Candidate method presumptive positive outcomes divided by the total number of trials

⁷POD_{CC} = Candidate method confirmed positive outcomes divided by the total number of trials

⁸dPOD_{CP} = Difference between the candidate method presumptive result and candidate method confirmed result POD values

⁹95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level

Table 39. Ceramic Tile Environmental Swabs Results: Analyzed on the Agilent AriaMx by the BACGene Evaluation Sheet¹, Candidate vs. Reference – POD Results (1)

STEC Strain	PREraser Treatment	Enrichment Time Point	CFU/ Test Area	N ^b	Candidate ²			Reference			dPOD _c ^f	95% CI ^g
					x ^c	POD _c ^d	95% CI	x ^c	POD _R ^e	95% CI		
<i>E. coli</i> O157:H7 ATCC ³ 43895	Without PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	10	0.50	0.30, 0.70	8	0.40	0.22, 0.61	0.10	-0.19, 0.37
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	11	0.55	0.34, 0.74	8	0.40	0.22, 0.61	0.15	-0.15, 0.41
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	10	0.50	0.30, 0.70	8	0.40	0.22, 0.61	0.10	-0.19, 0.37
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	12	0.60	0.39, 0.78	8	0.40	0.22, 0.61	0.20	-0.10, 0.46
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹The BACGene Evaluation Sheet and FastFinder software produced the same results.²The BACGene STEC workflow and the BACGene *E. coli* O157:H7 workflow produced the same results.³ATCC – American Type Culture Collection⁴MPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval⁵N = Number of test portions⁶x = Number of positive test portions⁷POD_c = Candidate method confirmed positive outcomes divided by the total number of trials⁸POD_R = Reference method confirmed positive outcomes divided by the total number of trials⁹dPOD_c = Difference between the confirmed candidate method result and reference method confirmed result POD values¹⁰95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level**Table 40. Ceramic Tile Environmental Swabs Results: Analyzed on the Bio-Rad CFX96 Touch Deep Well by the BACGene Evaluation Sheet, Candidate vs. Reference – POD Results (1)**

STEC Strain	PREraser Treatment	Enrichment Time Point	CFU/ Test Area	N ^b	Candidate ¹			Reference			dPOD _c ^f	95% CI ^g
					x ^c	POD _c ^d	95% CI	x ^c	POD _R ^e	95% CI		
<i>E. coli</i> O157:H7 ATCC ² 43895	Without PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	10	0.50	0.30, 0.70	8	0.40	0.22, 0.61	0.10	-0.19, 0.37
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	12	0.60	0.39, 0.78	8	0.40	0.22, 0.61	0.20	-0.10, 0.46
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	10	0.50	0.30, 0.70	8	0.40	0.22, 0.61	0.10	-0.19, 0.37
			570	5	4	0.80	0.38, 1.00	5	1.00	0.57, 1.00	-0.20	-0.62, 0.28
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	12	0.60	0.39, 0.78	8	0.40	0.22, 0.61	0.20	-0.10, 0.46
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹The BACGene STEC workflow and the BACGene *E. coli* O157:H7 workflow produced the same results.²ATCC – American Type Culture Collection³MPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval⁴N = Number of test portions⁵x = Number of positive test portions⁶POD_c = Candidate method confirmed positive outcomes divided by the total number of trials⁷POD_R = Reference method confirmed positive outcomes divided by the total number of trials⁸dPOD_c = Difference between the confirmed candidate method result and reference method confirmed result POD values⁹95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level

Table 41. Ceramic Tile Environmental Swabs Results: Analyzed on the Bio-Rad CFX96 Touch Deep Well by the FastFinder software, Candidate vs. Reference – POD Results (1)

STEC Strain	PREraser Treatment	Enrichment Time Point	CFU/ Test Area	N ^b	Candidate ¹			Reference			dPOD _c ^f	95% CI ^g
					x ^c	POD _c ^d	95% CI	x ^c	POD _R ^e	95% CI		
<i>E. coli</i> O157:H7 ATCC ² 43895	Without PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	10	0.50	0.30, 0.70	8	0.40	0.22, 0.61	0.10	-0.19, 0.37
			570	5	4	0.80	0.38, 1.00	5	1.00	0.57, 1.00	-0.20	-0.62, 0.28
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	12	0.60	0.39, 0.78	8	0.40	0.22, 0.61	0.20	-0.10, 0.46
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	9	0.45	0.26, 0.66	8	0.40	0.22, 0.61	0.05	-0.24, 0.33
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	12	0.60	0.39, 0.78	8	0.40	0.22, 0.61	0.20	-0.10, 0.46
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹The BACGene STEC workflow and the BACGene *E. coli* O157:H7 workflow produced the same results.²ATCC – American Type Culture Collection³MPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval⁴N = Number of test portions⁵x = Number of positive test portions⁶POD_c = Candidate method confirmed positive outcomes divided by the total number of trials⁷POD_R = Reference method confirmed positive outcomes divided by the total number of trials⁸dPOD_c = Difference between the confirmed candidate method result and reference method confirmed result POD values⁹95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level**Table 42. Ceramic Tile Environmental Swabs Results: Analyzed on the Agilent AriaMx by the BACGene Evaluation Sheet, Presumptive vs. Confirmed – POD Results (1)**

STEC Strain	PREraser Treatment	Enrichment Time Point	CFU/ Test Area	N ^b	Presumptive ¹			Confirmed ²			dPOD _{CP} ^f	95% CI ^g
					x ^c	POD _{CP} ^d	95% CI	x ^c	POD _{CC} ^e	95% CI		
<i>E. coli</i> O157:H7 ATCC ³ 43895	Without PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	12	0.60	0.39, 0.78	12	0.60	0.39, 0.78	0.00	-0.13, 0.13
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	12	0.60	0.39, 0.78	12	0.60	0.39, 0.78	0.00	N/A
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	10	0.50	0.30, 0.70	12	0.60	0.39, 0.78	-0.10	-0.28, 0.08
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	12	0.60	0.39, 0.78	12	0.60	0.39, 0.78	0.00	N/A
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹The BACGene STEC workflow and the BACGene *E. coli* O157:H7 workflow produced the same results.²Alternative and traditional confirmation procedures produced the same results.³ATCC – American Type Culture Collection⁴MPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval⁵N = Number of test portions⁶x = Number of positive test portions⁷POD_{CP} = Candidate method presumptive positive outcomes divided by the total number of trials⁸POD_{CC} = Candidate method confirmed positive outcomes divided by the total number of trials⁹dPOD_{CP} = Difference between the candidate method presumptive result and candidate method confirmed result POD values¹⁰95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level

Table 43. Ceramic Tile Environmental Swabs Results: Analyzed on the Agilent AriaMx by the FastFinder software, Presumptive vs. Confirmed – POD Results (1)

STEC Strain	PREraser Treatment	Enrichment Time Point	CFU/ Test Area	N ^b	Presumptive ¹			Confirmed ²			dPOD _{CP} ^f	95% CI ^g
					x ^c	POD _{CP} ^d	95% CI	x ^c	POD _{CC} ^e	95% CI		
<i>E. coli</i> O157:H7 ATCC ³ 43895	Without PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	-0.47, 0.47
			60	20	12	0.60	0.39, 0.78	12	0.60	0.39, 0.78	0.00	-0.13, 0.13
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	-0.47, 0.47
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	-0.47, 0.47
			60	20	13	0.65	0.43, 0.82	12	0.60	0.39, 0.78	0.05	-0.11, 0.21
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	-0.47, 0.47
	With PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	-0.47, 0.47
			60	20	10	0.50	0.30, 0.70	12	0.60	0.39, 0.78	-0.10	-0.28, 0.08
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	-0.47, 0.47
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	-0.47, 0.47
			60	20	12	0.60	0.39, 0.78	12	0.60	0.39, 0.78	0.00	-0.13, 0.13
			570	5	4	1.00	0.80, 1.00	5	1.00	0.57, 1.00	-0.20	-0.76, 0.36

¹The BACGene STEC workflow and the BACGene *E. coli* O157:H7 workflow produced the same results.²Alternative and traditional confirmation procedures produced the same results.³ATCC – American Type Culture Collection⁴MPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval⁵N = Number of test portions⁶x = Number of positive test portions⁷POD_{CP} = Candidate method presumptive positive outcomes divided by the total number of trials⁸POD_{CC} = Candidate method confirmed positive outcomes divided by the total number of trials⁹dPOD_{CP} = Difference between the candidate method presumptive result and candidate method confirmed result POD values¹⁰95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level**Table 44. Ceramic Tile Environmental Swabs Results: Analyzed on the Bio-Rad CFX96 Touch Deep Well by the BACGene Evaluation Sheet, Presumptive vs. Confirmed – POD Results (1)**

STEC Strain	PREraser Treatment	Enrichment Time Point	CFU/ Test Area	N ^b	Presumptive ¹			Confirmed ²			dPOD _{CP} ^f	95% CI ^g
					x ^c	POD _{CP} ^d	95% CI	x ^c	POD _{CC} ^e	95% CI		
<i>E. coli</i> O157:H7 ATCC ³ 43895	Without PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	12	0.60	0.39, 0.78	12	0.60	0.39, 0.78	0.00	N/A
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	13	0.65	0.43, 0.82	12	0.60	0.39, 0.78	0.05	-0.11, 0.21
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	10	0.50	0.30, 0.70	12	0.60	0.39, 0.78	-0.10	-0.28, 0.08
			570	5	4	1.00	0.80, 1.00	5	1.00	0.57, 1.00	-0.20	-0.76, 0.36
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	12	0.60	0.39, 0.78	12	0.60	0.39, 0.78	0.00	N/A
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹The BACGene STEC workflow and the BACGene *E. coli* O157:H7 workflow produced the same results.²Alternative and traditional confirmation procedures produced the same results.³ATCC – American Type Culture Collection⁴MPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval⁵N = Number of test portions⁶x = Number of positive test portions⁷POD_{CP} = Candidate method presumptive positive outcomes divided by the total number of trials⁸POD_{CC} = Candidate method confirmed positive outcomes divided by the total number of trials⁹dPOD_{CP} = Difference between the candidate method presumptive result and candidate method confirmed result POD values¹⁰95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level

Table 45. Ceramic Tile Environmental Swabs Results: Analyzed on the Bio-Rad CFX96 Touch™ Deep Well by the Fast Finder software, Presumptive vs. Confirmed – POD Results (1)

STEC Strain	PREraser Treatment	Enrichment Time Point	CFU/ Test Area	N ^b	Presumptive ¹			Confirmed ²			dPOD _{CP} ^f	95% CI ^g
					x ^c	POD _{CP} ^d	95% CI	x ^c	POD _{CC} ^e	95% CI		
<i>E. coli</i> O157:H7 ATCC ³ 43895	Without PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	12	0.60	0.39, 0.78	12	0.60	0.39, 0.78	0.00	N/A
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	13	0.65	0.43, 0.82	12	0.60	0.39, 0.78	0.05	-0.11, 0.21
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	10	0.50	0.30, 0.70	12	0.60	0.39, 0.78	-0.10	-0.28, 0.08
			570	5	4	1.00	0.80, 1.00	5	1.00	0.57, 1.00	-0.20	-0.76, 0.36
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	12	0.60	0.39, 0.78	12	0.60	0.39, 0.78	0.00	N/A
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹The BACGene STEC workflow and the BACGene *E. coli* O157:H7 workflow produced the same results.²Alternative and traditional confirmation procedures produced the same results.³ATCC – American Type Culture Collection^aMPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval^bN = Number of test portions^cx = Number of positive test portions^dPOD_{CP} = Candidate method presumptive positive outcomes divided by the total number of trials^ePOD_{CC} = Candidate method confirmed positive outcomes divided by the total number of trials^fdPOD_{CP} = Difference between the candidate method presumptive result and candidate method confirmed result POD values^g95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level**Table 46. Ceramic Tile Environmental Swabs Results: Analyzed on the Agilent AriaMx™ by the BACGene Evaluation Sheet, Candidate vs. Reference – POD Results (1)**

STEC Strain	PREraser Treatment	Enrichment Time Point	CFU/ Test Area	N ^b	Candidate			Reference			dPOD _C ^f	95% CI ^g
					x ^c	POD _C ^d	95% CI	x ^c	POD _R ^e	95% CI		
<i>E. coli</i> O145 PSU ¹ 7.1711	Without PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	8	0.40	0.22, 0.61	6	0.30	0.15, 0.52	0.10	-0.18, 0.36
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	8	0.40	0.22, 0.61	6	0.30	0.15, 0.52	0.10	-0.18, 0.36
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	8	0.40	0.22, 0.61	6	0.30	0.15, 0.52	0.10	-0.18, 0.36
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	8	0.40	0.22, 0.61	6	0.30	0.15, 0.52	0.10	-0.18, 0.36
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹PSU – Pennsylvania State University Culture Collection^aMPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval^bN = Number of test portions^cx = Number of positive test portions^dPOD_C = Candidate method confirmed positive outcomes divided by the total number of trials^ePOD_R = Reference method confirmed positive outcomes divided by the total number of trials^fdPOD_C = Difference between the confirmed candidate method result and reference method confirmed result POD values^g95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level

Table 47. Ceramic Tile Environmental Swabs Results: Analyzed on the Agilent AriaMx™ by the FastFinder software, Candidate vs. Reference – POD Results (1)

STEC Strain	PREraser Treatment	Enrichment Time Point	CFU/ Test Area	N ^b	Candidate			Reference			dPOD _c ^f	95% CI ^g
					x ^c	POD _c ^d	95% CI	x ^c	POD _R ^e	95% CI		
<i>E. coli</i> O145 PSU ¹ 7.1711	Without PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	8	0.40	0.22, 0.61	6	0.30	0.15, 0.52	0.10	-0.18, 0.36
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	8	0.40	0.22, 0.61	6	0.30	0.15, 0.52	0.10	-0.18, 0.36
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	8	0.40	0.22, 0.61	6	0.30	0.15, 0.52	0.10	-0.18, 0.36
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	7	0.35	0.18, 0.57	6	0.30	0.15, 0.52	0.05	-0.23, 0.32
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹PSU – Pennsylvania State University Culture Collection^aMPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval^bN = Number of test portions^cx = Number of positive test portions^dPOD_c = Candidate method confirmed positive outcomes divided by the total number of trials^ePOD_R = Reference method confirmed positive outcomes divided by the total number of trials^fdPOD_c = Difference between the confirmed candidate method result and reference method confirmed result POD values^g95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level**Table 48. Ceramic Tile Environmental Swabs Results: Analyzed on the Bio-Rad CFX96 Touch™ Deep Well by the BACGene Evaluation Sheet, Candidate vs. Reference – POD Results (1)**

STEC Strain	PREraser Treatment	Enrichment Time Point	CFU/ Test Area	N ^b	Candidate			Reference			dPOD _c ^f	95% CI ^g
					x ^c	POD _c ^d	95% CI	x ^c	POD _R ^e	95% CI		
<i>E. coli</i> O145 PSU ¹ 7.1711	Without PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	8	0.40	0.22, 0.61	6	0.30	0.15, 0.52	0.10	-0.18, 0.36
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	8	0.40	0.22, 0.61	6	0.30	0.15, 0.52	0.10	-0.18, 0.36
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	8	0.40	0.22, 0.61	6	0.30	0.15, 0.52	0.10	-0.18, 0.36
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	8	0.40	0.22, 0.61	6	0.30	0.15, 0.52	0.10	-0.18, 0.36
			570	5	4	0.80	0.38, 1.00	5	1.00	0.57, 1.00	-0.20	-0.62, 0.28

¹PSU – Pennsylvania State University Culture Collection^aMPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval^bN = Number of test portions^cx = Number of positive test portions^dPOD_c = Candidate method confirmed positive outcomes divided by the total number of trials^ePOD_R = Reference method confirmed positive outcomes divided by the total number of trials^fdPOD_c = Difference between the confirmed candidate method result and reference method confirmed result POD values^g95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level

Table 49. Ceramic Tile Environmental Swabs Results: Analyzed on the Bio-Rad CFX96 Touch™ Deep Well by the FastFinder software, Candidate vs. Reference – POD Results (1)

STEC Strain	PREraser Treatment	Enrichment Time Point	CFU/ Test Area	N ^b	Candidate			Reference			dPOD _c ^f	95% CI ^g
					x ^c	POD _c ^d	95% CI	x ^c	POD _R ^e	95% CI		
<i>E. coli</i> O145 PSU ¹ 7.1711	Without PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	8	0.40	0.22, 0.61	6	0.30	0.15, 0.52	0.10	-0.18, 0.36
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	8	0.40	0.22, 0.61	6	0.30	0.15, 0.52	0.10	-0.18, 0.36
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	7	0.35	0.18, 0.57	6	0.30	0.15, 0.52	0.05	-0.23, 0.32
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	8	0.40	0.22, 0.61	6	0.30	0.15, 0.52	0.10	-0.18, 0.36
			570	5	4	0.80	0.38, 1.00	5	1.00	0.57, 1.00	-0.20	-0.62, 0.28

¹PSU – Pennsylvania State University Culture Collection^aMPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval^bN = Number of test portions^cx = Number of positive test portions^dPOD_c = Candidate method confirmed positive outcomes divided by the total number of trials^ePOD_R = Reference method confirmed positive outcomes divided by the total number of trials^fdPOD_c = Difference between the confirmed candidate method result and reference method confirmed result POD values^g95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level**Table 50. Ceramic Tile Environmental Swabs Results: Analyzed on the Agilent AriaMx™ by the BACGene Evaluation Sheet, Presumptive vs. Confirmed – POD Results (1)**

STEC Strain	PREraser Treatment	Enrichment Time Point	CFU/ Test Area	N ^b	Presumptive			Confirmed ¹			dPOD _{CP} ^f	95% CI ^g
					x ^c	POD _{CP} ^d	95% CI	x ^c	POD _{CC} ^e	95% CI		
<i>E. coli</i> O145 PSU ² 7.1711	Without PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	9	0.45	0.26, 0.66	8	0.40	0.22, 0.61	0.05	-0.11, 0.21
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	12	0.60	0.39, 0.78	8	0.40	0.22, 0.61	0.20	-0.01, 0.41
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	8	0.40	0.22, 0.61	8	0.40	0.22, 0.61	0.00	N/A
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	8	0.40	0.22, 0.61	8	0.40	0.22, 0.61	0.00	N/A
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹Alternative and traditional confirmation procedures produced the same results.²PSU – Pennsylvania State University Culture Collection^aMPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval^bN = Number of test portions^cx = Number of positive test portions^dPOD_{CP} = Candidate method presumptive positive outcomes divided by the total number of trials^ePOD_{CC} = Candidate method confirmed positive outcomes divided by the total number of trials^fdPOD_{CP} = Difference between the candidate method presumptive result and candidate method confirmed result POD values^g95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level

Table 51. Ceramic Tile Environmental Swabs Results: Analyzed on the Agilent AriaMx™ by the FastFinder software, Presumptive vs. Confirmed – POD Results (1)

STEC Strain	PREraser Treatment	Enrichment Time Point	CFU/ Test Area	N ^b	Presumptive			Confirmed ¹			dPOD _{CP} ^f	95% CI ^g
					x ^c	POD _{CP} ^d	95% CI	x ^c	POD _{CC} ^e	95% CI		
<i>E. coli</i> O145 PSU ² 7.1711	Without PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	9	0.45	0.26, 0.66	8	0.40	0.22, 0.61	0.05	-0.11, 0.21
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	12	0.60	0.39, 0.78	8	0.40	0.22, 0.61	0.20	-0.01, 0.41
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	8	0.40	0.22, 0.61	8	0.40	0.22, 0.61	0.00	N/A
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	7	0.35	0.18, 0.57	8	0.40	0.22, 0.61	-0.05	-0.21, 0.11
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A

¹Alternative and traditional confirmation procedures produced the same results.²PSU – Pennsylvania State University Culture Collection³MPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval⁴N = Number of test portions⁵x = Number of positive test portions⁶POD_{CP} = Candidate method presumptive positive outcomes divided by the total number of trials⁷POD_{CC} = Candidate method confirmed positive outcomes divided by the total number of trials⁸dPOD_{CP} = Difference between the candidate method presumptive result and candidate method confirmed result POD values⁹95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level**Table 52. Ceramic Tile Environmental Swabs Results: Analyzed on the Bio-Rad CFX96 Touch™ Deep Well by the BACGene Evaluation Sheet, Presumptive vs. Confirmed – POD Results (1)**

STEC Strain	PREraser Treatment	Enrichment Time Point	CFU/ Test Area	N ^b	Presumptive			Confirmed ¹			dPOD _{CP} ^f	95% CI ^g
					x ^c	POD _{CP} ^d	95% CI	x ^c	POD _{CC} ^e	95% CI		
<i>E. coli</i> O145 PSU ² 7.1711	Without PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	9	0.45	0.26, 0.66	8	0.40	0.22, 0.61	0.05	-0.11, 0.21
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	12	0.60	0.39, 0.78	8	0.40	0.22, 0.61	0.20	-0.01, 0.41
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	8	0.40	0.22, 0.61	8	0.40	0.22, 0.61	0.00	N/A
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	8	0.40	0.22, 0.61	8	0.40	0.22, 0.61	0.00	N/A
			570	5	4	0.80	0.38, 1.00	5	1.00	0.57, 1.00	-0.20	-0.76, 0.36

¹Alternative and traditional confirmation procedures produced the same results.²PSU – Pennsylvania State University Culture Collection³MPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval⁴N = Number of test portions⁵x = Number of positive test portions⁶POD_{CP} = Candidate method presumptive positive outcomes divided by the total number of trials⁷POD_{CC} = Candidate method confirmed positive outcomes divided by the total number of trials⁸dPOD_{CP} = Difference between the candidate method presumptive result and candidate method confirmed result POD values⁹95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level

Table 53. Ceramic Tile Environmental Swabs Results: Analyzed on the Bio-Rad CFX96 Touch™ Deep Well by the FastFinder software, Presumptive vs. Confirmed – POD Results (1)

STEC Strain	PREraser Treatment	Enrichment Time Point	CFU/ Test Area	N ^b	Presumptive			Confirmed ¹			dPOD _{CP} ^f	95% CI ^g
					x ^c	POD _{CP} ^d	95% CI	x ^c	POD _{CC} ^e	95% CI		
<i>E. coli</i> O145 PSU ² 7.1711	Without PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	10	0.50	0.30, 0.70	8	0.40	0.22, 0.61	0.10	-0.08, 0.28
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	12	0.60	0.39, 0.78	8	0.40	0.22, 0.61	0.20	-0.01, 0.41
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
	With PREraser	10 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	7	0.35	0.18, 0.57	8	0.40	0.22, 0.61	-0.05	-0.21, 0.11
			570	5	5	1.00	0.57, 1.00	5	1.00	0.57, 1.00	0.00	N/A
		24 Hour	-	5	0	0.00	0.00, 0.43	0	0.00	0.00, 0.43	0.00	N/A
			60	20	8	0.40	0.22, 0.61	8	0.40	0.22, 0.61	0.00	N/A
			570	5	4	0.80	0.38, 1.00	5	1.00	0.57, 1.00	-0.20	-0.76, 0.36

¹Alternative and traditional confirmation procedures produced the same results.²PSU – Pennsylvania State University Culture Collection³MPN = Most Probable Number is calculated using the LCF MPN calculator ver. 1.6 provided by AOAC RI, with 95% confidence interval⁴N = Number of test portions⁵x = Number of positive test portions⁶POD_{CP} = Candidate method presumptive positive outcomes divided by the total number of trials⁷POD_{CC} = Candidate method confirmed positive outcomes divided by the total number of trials⁸dPOD_{CP} = Difference between the candidate method presumptive result and candidate method confirmed result POD values⁹95% CI = If the confidence interval of a dPOD does not contain zero, then the difference is statistically significant at the 5% level

DISCUSSION OF MODIFICATION APPROVED JANUARY 2021 (5)

For the BACGene real-time PCR kits for the detection of bacterial pathogens in food and feed, as well as in environmental samples, the two tested alternative non-ionic detergents, Ecosurf™ SA-9 and Tween® 20, show a very similar lysis efficiency in comparison to the reference detergent Triton® X-100.

The product consistency study showed good agreement between the three lysis buffer formulations each tested with three different lots, with consistent and reproducible results. Cq-values in spiked matrixes were constant and reproducible, with mean Cq-values CV < 5% and thus it does not indicate a systematic inconsistency from production of individual lots with the different detergents. For BACGene *Salmonella* spp. 2, two inhibited samples were observed for Lot 2 and one inhibited sample for Lot 3. For lysis efficiency comparison those three replicates were taken out. Since all other replicates for those two lots are fine, most probably wrong handling (taking the lysate not from the upper half) can be seen as the reason for that.

Additional to the product consistency study, a matrix test was performed with the three different lysis buffer types, each with three relevant food and environmental samples, which were inoculated with the respective target strains. Enrichment was performed according to the instructions and lysis was performed with the certified buffer containing Triton® X-100 in comparison to the two alternatives, Ecosurf™ SA-9 and Tween® 20.

In none of the tested samples for none of the kits there was a significant difference between the lysis efficiencies observed. Cq-values are very similar between the reference condition and both alternatives.

Therefore, it may be concluded that this study demonstrates a similar performance of both tested alternatives, Ecosurf™ SA-9 and Tween® 20, compared to Triton® X-100.

DISCUSSION OF MODIFICATION APPROVED JANUARY 2024 (6)

In total, 2457, 5103 and 3276 datasets have been re-evaluated for the BACGene *Salmonella* spp. method, the three BACGene *Listeria* methods and the BACGene *E. coli* O157:H7 and STEC TOP7 workflows, respectively. Within the tested inclusivity/exclusivity datasets, no deviation to the initial validation data was observed. These results strongly indicate that there is no negative impact of the usage of PURE on the specificity of any tested method.

In matrix datasets, a total of 7 deviations regarding final results occurred: 3 for the BACGene *Salmonella* spp. method, 1 for the three BACGene *Listeria* methods and 3 for the BACGene *E. coli* O157:H7 and the STEC TOP7 workflows. In the deviating samples, a fluorescence increase at some point during the PCR reaction was visible in the raw data. The change in final scoring of these samples from presumptive negative to presumptive positive, was subsequently due to a different interpretation/processing of the respective data by either PURE or the real-time PCR cyclor software, which is the basis for the Excel-based data evaluation. For such borderline signals, it is not unexpected that different algorithms yield different results in data interpretation.

Most important, the observed deviations in single sample results, did not have a relevant negative impact on the overall performance characteristics of the evaluated methods.

REFERENCES CITED

1. Klass, N., Bastin, B., Benzinger, Jr., M.B., Agin, J., Bleichner, L., and Bahrtdt, C., Validation of the BACGene, *E. coli* O157:H7 Workflow and BACGene STEC TOP7 Workflow, AOAC Performance Tested MethodsSM certification number 022002.
2. U. S. Department of Agriculture-Food Safety and Inspection Service Microbiology Laboratory Guidebook (MLG), 5C.00, *Detection, Isolation and Identification of Top Seven Shiga Toxin-Producing Escherichia coli (STECs) from Meat Products and Carcass and Environmental Sponges*. (Accessed November 2019) <https://www.fsis.usda.gov/wps/wcm/connect/7ffc02b5-3d33-4a79-b50c-81f208893204/MLG-5B.pdf?MOD=AJPERES> (2)
3. Food and Drug Administration (FDA) Bacteriological Analytical Manual (BAM), Chapter 4A, *Diarrheagenic Escherichia coli*. (Accessed December 2019) <https://www.fda.gov/food/laboratory-methods-food/bam-diarrheagenic-escherichia-coli>
4. Varga, S. and Horváth, N., Validation of the Additional Manufacturing Location for BACGene *Listeria* spp., AOAC Performance Tested MethodsSM certification number 022002. Approved November 2020.
5. Bleichner, L. and Bleicher, V., Validation of the replacement of Triton X-100 in the BACGene Lysis Buffer, AOAC Performance Tested MethodsSM certification number 022002. Approved January 2021.
6. Bleichner, L. and Kohl, S., Validation of Gold Standard Diagnostics Freiburg GmbH BACGene Test Kits for the Extension for PURE evaluation software BACGene *Salmonella* spp. AOAC Performance Tested MethodSM 121501, BACGene *Listeria* Multiplex AOAC Performance Tested MethodSM 061701, BACGene *Listeria* spp. AOAC Performance Tested MethodSM 061702, BACGene *Listeria monocytogenes* AOAC Performance Tested MethodSM 061703, BACGene *E. coli* O157:H7 workflow AOAC Performance Tested MethodSM 022002, BACGene STEC TOP7 workflow AOAC Performance Tested MethodSM 022003, AOAC Performance Tested MethodsSM certification number 022002. Approved January 2024.